



CHILLED WATER CIRCULATOR - 4L

Applications:

- ♦ Ideally suitable for chilling and condensing operations. Ideal companion for Rotary Vacuum Evaporator, Viscosity, Spectrophotometry, Water Baths upto 28 lit and for general lab requirements of chilling in closed loop systems.

EQUITRON Advantages:

- ♦ Environment friendly CFC free hermetically sealed compressor.
- ♦ Sleek & compact design with easy maneuverability.
- ♦ Can be placed under workbench.

Construction / Standard Features:

- ♦ Stainless steel base & top. Panels of steel - powder coated.
- ♦ Four castors for easy movement.
- ♦ Supplied with 1 met x 2 Silicone tubings (8x11 mm) with insulation.
- ♦ In-built Digital Temperature Controller cum Indicator.
- ♦ Cooling capacity 750 W at 15°C; min. range upto 0°C (at ambient 25°C).
- ♦ Reservoir capacity of approximately 4 lit.
- ♦ Circulation pump of app 15 LPM at zero head - suitable for closed loop systems.



External size (L/W/D): 56 x 41 x 64 cm

Model	Description	Cooling Capacity	App. Shipping Details Net / Gr. / Volume
#8502.ROT.000	Chilled Water Circulator - 4L	750 W at 15°C	41 kg / 65 kg / 0.345 m ³

CHILLED WATER CIRCULATOR - 25L

Applications:

- ♦ Ideally suitable for chilling and condensing operations.
- ♦ Ideal companion for Bilster Pack Machines, 20 lit Rotary Vacuum Evaporators, and for general production requirements of chilling in closed loop systems.

EQUITRON Advantages:

- ♦ Environment friendly CFC free hermetically sealed rotary compressor.
- ♦ Sleek & compact design with easy maneuverability.

Construction / Standard Features:

- ♦ Sturdy angular welded frame - full stainless steel construction.
- ♦ Panels powder coated.
- ♦ Four heavy duty castors for easy movement.
- ♦ In-built Digital Temperature Controller cum Indicator.
- ♦ Temperature range upto +5°C (at ambient 25°C).
- ♦ Reservoir capacity of approximately 25 lit.
- ♦ Heavy duty water circulation pump of 180 W with 10 LPM flow rate at 15 met head.
- ♦ Suitable for closed loop systems.
- ♦ With low liquid level indication & automatic cut-off.



External size (L/W/D): 65 x 62 x 84 cm

Model	Description	Cooling Capacity	App. Shipping Details Net / Gr. / Volume
#8502.025.000	Chilled Water Circulator - 25L	2500 W at 15°C	81 kg / 110 kg / 0.58 m ³

MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013. Tamil Nadu, India
Ph: 044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz