



For zero infections....

EQUITRON®

CLASS-B TABLE TOP AUTOCLAVES

Applications:

- ♦ **EQUITRON** Table Top Vacuum High Speed Sterilizers are ideal for all applications requiring rapid, reliable and total destruction of all living microorganisms.
- ♦ Ideal for Ophthalmic Surgeons, Dental Surgeons, Surgical Doctors, Polyclinics and Hospitals, Operation theatres, QC / Micro as a stand-by unit.

EQUITRON Advantages:

- ♦ It adopts European Class B standard.
- ♦ With pulse pre-vacuum - essential for sterilizing all wrapped (single / double) or unwrapped items, including porous items that do not exceed the specifications of Hollow load Type A.
- ♦ A pre defined program list for easy operational settings: Reproducible sterilization cycles are accomplished through the defined programs.
- ♦ Electronic control module offers microprocessor controls and LED display of operating parameters, where temperature, pressure and sterilisation time are constantly displayed.
- ♦ A drying cycle helps prevent rusting of delicate instruments and helps dry clothing or absorbent materials.



Construction / Standard Features:

- ♦ Safer and more reliable door structure. A heat-insulated door with locking feature prevents opening while the chamber is under pressure to help minimize dangers to operator.
- ♦ Exhaust steam is discharged to an inbuilt water collection tank thus keeping the room steam free.

Technical Specifications:

Sterilizing Temperature	121°C or 134°C
Test	Bowie & Dick test
Safety Features	Pressure Release Valve and computerized diagnosing and alarm system
Water supply system	Reservoirs (2Ltr fresh water supply system and 1.6Ltr waste water reservoirs)
Voltage	220V±10% 50/60Hz, single phase
Power	2400W / 10A

Model	Product	Chamber Size	External Dimension W/D/H in mm	Heater (W)	Trays	App. Shipping Details (Air) Net / Gr. / Volume
#7471-B-22	EQUITRON "Class - B" Autoclave	Ø250mm x 450mm (22L approx.)	450 x 650+150 x 400	1500 Band heater 750 for Steam generator	3	45 kg / 72 kg / 0.35m ³

MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013, Tamil Nadu, India
Ph:044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz