

# COLONY COUNTER - GMP Model

### Applications:

- ◆ LAPIZ Digital Colony Counter is ideal for fast and fatigue free counting of Bacterial Colonies.

### LAPIZ Advantages:

- ◆ Circular and powerful magnifying lens - reduces operator fatigue.
- ◆ Streamlined design ensures fast and accurate count.
- ◆ Desk top compact model.

### Construction:

- ◆ Re-designed model with advanced sensors enables usage of any commercially available hard tipped marker pen.
- ◆ LED Display with audio beep at every count.
- ◆ Wouffhugel grid glass plate.
- ◆ On marking the colony by applying a slight pressure with the marker pen, a count is registered and confirmed by an audio beep.
- ◆ Full Stainless steel fabricated body duly heat cured epoxy coated from outside for better aesthetics.

### Features:

- ◆ Magnification: 2x.
- ◆ Maximum Count: 9,999.
- ◆ Manual count switch for counting without "Marker Pen".



Model	Product	Maximum Count	App. Shipping Details Net / Gr. / Volume
#0671M	Digital - Desk Model - SS GMP Model	9,999	4.5 kg / 5.0 kg (41 x 32 x 30 cm) 0.04 m <sup>3</sup>

### OPTIONAL ACCESSORIES FOR COLONY COUNTER

Model	Product	Details
#7960	Petri Platter	Manual for 100 mm plates
-#7963-	-Petri-Plate-Adaptors-	-Suitable for 50 and 75 mm Ø plates-

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.  
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical, Standard / PAD / Fully Automatic / Non Electric / Jacketed, Autoclave Portable, Standard / Automatic / Non Electric, Autoclave SLL Series, Vacuum, Autoclave SLE series, Table Top Autoclaves: Normal / Class B - Rotary Vacuum Evaporator - Vacuum Controller - Chilled Water Circulator - PTFE Diaphragm Vacuum Pump - BOD Incubator - BSC Clean Bench - Vacuum Oven - Water Still / Stone-It / Automatic Water Still - Water Bath: Sterile / Constant / (Serological) / Concentric Rings / Swaging / Refrigerated Circulating Bath - Round Bath: Clear Bottom / Immersion Series - Ultrasonic Bath: Digital / LCD - Dry Bath - Colony Counter - Hot plate with Magnetic Stirrer - Anaerobic Jar: SS / Polycarbonate - Pipette & Petri Dish Sterilizing Box - Test Tube Rack - ESR Westergren & Wintrobe Apparatus - Slide Staining Stand / Rack / Trough - Slide Trays - Slide Cabinet

## MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,  
Royapuram, Chennai - 600 013. Tamil Nadu, India  
Ph:044-25984175, 45526026, 094441 62346  
Email: multilabindia@gmail.com, www.multilab.biz