



REFRIGERATED CIRCULATING BATH 13 & 22 Ltr.

Applications:

- ♦ Ideal for Rotational Viscosity tests, Spectrophotometers, Refractometers, Electrophoresis, Isoelectric Focusing, Calibration, Chromatography, Condensers, Densitometers, Polarimeters and General Laboratory Cooling.

EQUIBATH advantages:

- ♦ EQUIBATH provides heating & cooling.
- ♦ Eliminates all the unwieldiness and hassles of using a separate Bath & Chiller.
- ♦ Built in circulator pumps to external devices.
- ♦ It offers highest level of performance, flexibility and control for the most demanding applications.
- ♦ Quiet operation.
- ♦ Wide operating range.

Construction:

- ♦ Corrosion resistant stainless steel bath & lid.
- ♦ Exterior construction of stainless steel with heat cured epoxy coating.

Standard Features:

- ♦ Temperature Range -10 to +120°C.
- ♦ PID Temperature controller with timer, alarm & over temperature protection.
- ♦ Control accuracy of $\pm 0.1^\circ\text{C}$ at -10°C .
- ♦ Internal circulation enables precise temperature uniformity in the bath better than $\pm 0.3^\circ\text{C}$.
- ♦ Resolution 0.1°C .
- ♦ Impeller type external circulation pump with 8 Lpm capacity.
- ♦ Cooling capacity 350W at 20°C .



Model	Bath	Temp. Range	Capacity	Heater	Usable Space L/W/H in mm	App. Shipping Details Net / Gr. / Volume
#8508.13L	Refrigerated Circulating EQUIBATH	-10 to +120°C	13 lit.	750 W	240 x 150 x 150	35 kg / 45 kg / 0.15m ³
#8508.22L	Refrigerated Circulating EQUIBATH	-10 to +120°C	22 lit.	1300 W	320 x 250 x 150	58 kg / 75 kg / 0.25m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PID / Fully Automatic / Jacketed, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series - Media Preparator - Table Top Autoclaves: Normal / Class B - UV Sterilization Cabinet - Rotary Vacuum Evaporator - Vacuum Controller - Chilled Water Circulator - PFE Diaphragm Vacuum Pump - BOD Incubator - BSC Clean Bench - Water Still / Stone-In / Automatic Water Still - Water Bath: Stirred (Constant) / Serological / Concentric Rings / Shaking / Refrigerated Circulating Bath - Round Bath: Clear Bottom / Immersion Series / ECO Series - Ultrasonic Bath: Digital / LCD - Dry Bath - Colony Counter - Hot Plate with Magnetic Stirrer - Anaerobic Jar: SS / Polycarbonate - Pipette & Petri Dish Sterilizing Box - Test Tube Rack - ESR Westergren & Wintrobe Apparatus - Slide Staining Stand / Rack / Trough - Slide Trays - Slide Cabinet

MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013, Tamil Nadu, India
Ph:044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz