



For convenience and economy

**EQUITRON**<sup>®</sup>

## **ROUND WATER / OIL BATH**

### **IMMERSION SERIES**

#### **Digital**

#### **Applications:**

- ❖ **EQUITRON** Round Baths - with Digital control are ideal for reactions, extractions, analysis and tests.

#### **EQUITRON Advantages:**

- ❖ Double walled construction with efficient insulation system.
- ❖ Bath chamber of stainless steel.
- ❖ Outer body and control panel of MS duly finished in heat cured epoxy coating.
- ❖ On/off switch for mains and heater indication - on the front panel.
- ❖ Stainless Steel immersion heater - peripheral shape to eliminate direct contact with flask - even if flask touches the bottom.



#### **Choice of Control:**

- ❖ Digital Range: Ambient +10 to 250°C.
- ❖ Control accuracy: ±2°C or better.

Model	Round Flask Capacity	Bath size Ø x deep in mm	Temperature Range Water cum Oil - DWO	Heater	App. Shipping Details Net / Gr. / Volume
#8481SDWO	1 lit.	Ø 165 x 90	Ambient + 10°C to 250°C	400 W	3.5 kg / 4.0 kg / 0.020 m <sup>3</sup>
#8483SDWO	3 lit.	Ø 230 x 130		750 W	5.0 kg / 6.0 kg / 0.025 m <sup>3</sup>
#8485SDWO	5 lit.	Ø 260 x 150		1000 W	7.0 kg / 8.0 kg / 0.030 m <sup>3</sup>
#8490SDWO	10 lit.	Ø 350 x 200		1300 W	15.0 kg / 17.0 kg / 0.06 m <sup>3</sup>

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.  
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical / Standard / PAD / Fully Automatic / Jacketed, Autoclave Portable / Standard / Automatic, Autoclave SLL Series, Vacuum, Autoclave SLE series - Media Preparator / Table Top Autoclaves: Normal / Class B - UV Sterilization Cabinet / Rotary Vacuum Evaporator / Vacuum Controller / Chilled Water Circulator / PTFE Diaphragm Vacuum Pump / BOD Incubator / BSC Clean Bench - Water Still / Stone-free / Automatic Water Still - Water Bath: Stirred / Constant / Serological / Concentric Rings / Shaking / Refrigerated Circulating Bath - Round Bath: Clear Bottom / Immersion Series / ECO Series - Ultrasonic Bath: Digital / LCD - Dry Bath - Counting Counter - Hot Plate with Magnetic Stirrer - Anaerobic Jar: SS / Polycarbonate - Pipette & Petri Dish Sterilizing Box - Test Tube Rack - BSK Westergren & Wintrabe Apparatus - Slide Staining Stand/Rack / Trough / Slide Trays - Slide Cabinet

**MultiLab**

No.51/23, Arathoon Road, Off Appavu Lane,  
Royapuram, Chennai - 600 013, Tamil Nadu, India  
Ph:044-25984175, 45526026, 094441 62346  
Email: multilabindia@gmail.com, www.multilab.biz