



SEROLOGICAL WATER BATH - THERMOSTATIC

EQUITRON Advantages:

- ◆ Fully stainless steel, including bottom, makes it rust proof.
- ◆ Easy to set thermostat with 37°C calibrated mark.
- ◆ Ready to use - supplied with thermometer.

Construction:

- ◆ Double walled full stainless steel construction including bottom.
- ◆ Glasswool insulation for energy efficiency.

Standard Features:

- ◆ Capillary thermostat with $\pm 0.5^\circ\text{C}$ accuracy at 37°C.
- ◆ Energy efficient tubular heater.
- ◆ Range: Ambient plus 5°C to 65°C.
- ◆ Supplied with Thermometer, Clip and Mains Cord.



Model	Size L/W/D in cm	Racks (optional) Kahn 30 #6295	Kahn 48 #6299	Optional SS Gabled Lid Model	Optional Float Ball Packs Model	Heater	App. Shipping Details Net / Gr. / Volume
#6701TH	27 x 13 x 12	-	1	#6901	#6990-50 x 1	500 W	3.5 kg / 4.0 kg / 0.030 m ³
#6721TH	30 x 27 x 12	1	3	#6921	#6990-50 x 2	500 W	5.5 kg / 6.1 kg / 0.050 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

ACRYLIC WATER BATH 5L - THERMOSTATIC

EQUITRON Advantages:


- ◆ Acrylic body for full visibility of samples.
- ◆ Totally rust proof construction.

Construction:

- ◆ Double jointed Acrylic body for leak proof construction.
- ◆ SS rust proof control panel construction.

Standard Features:

- ◆ Temperature range ambient +5°C to 60°C.
- ◆ Supplied with thermometer & clip.
- ◆ Special design body for $\pm 0.5^\circ\text{C}$ control accuracy.
- ◆ Energy efficient top suspending tubular heater.



Model	Size L/W/D in cm	OPTIONAL Test Tube Rack - #6235-Acr-8	Heater	App. Shipping Details Net / Gr. / Volume
#6700TH	32 x 12 x 13	Ø20 x 8 Holes - 2 tier	500 W	2.8 kg / 3.1 kg / 0.037 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.