



For precise chemical & biochemical reactions

EQUITRON[®]

SHAKING BATH

Applications:

- ◆ **EQUITRON** Shaking Baths are ideal for all applications requiring precise and reproducible results for incubations, fermentations, homogenization, chemical and biochemical reactions, enzyme and tissue studies.

EQUITRON Advantages:

- ◆ Digital RPM as standard.
- ◆ Constant shaking speed even during heater operation.

Construction:

- ◆ Double walled, insulated, full stainless steel construction.
- ◆ Hinged Gabled SS Lid.
- ◆ Removable SS tray of modular design to accommodate flask holders of different capacities.
- ◆ Independent over-temperature self resetting safety cut-out with alarm.

Standard Features:

- ◆ Control: Digital 0.1°C resolution.
- ◆ Control Accuracy: $\pm 0.2^{\circ}\text{C}$ at 37°C.
- ◆ RPM: Digital.
- ◆ Temperature Range: Ambient + 5 to 80°C.
- ◆ Shaking Speed: 30 to 175 oscillations.
- ◆ Control Accuracy: ± 3 RPM.
- ◆ Stroke Length: 25 mm.



OPTIONAL ACCESSORIES

OPTIONAL ACCESSORIES		
#S3101-10	Flask holder for 100 ml flask	Max 9 nos.
#S3101-15	Flask holder for 500 ml flask	Max 4 nos.
#S3101-20	Flask holder for 1000 ml flask	Max 2 nos.
#S3101-25	Flask holder for 2000 ml flask	Max 1 no.
#S3101-50	Tray Universal for shaking bath	--

Model	Size L/W/D in cm	Flask holders supplied as Standard	Shaking Speed	App. Shipping Details Net / Gr. / Volume
#8406	45 x 28 x 21	6 nos. of 250 ml	30 to 175 osc.	32 kg / 49 kg / 0.27 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PID / Fully Automatic / Jacketed, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series - Media Preparator - Table Top Autoclaves: Normal / Class B - UV Sterilization Cabinet - Rotary Vacuum Evaporator - Vacuum Controller - Distill Water Circulator - FTFE Diaphragm Vacuum Pump - BOD Incubator - BSC Clean Bench - Water SSB / Steam-Box / Automatic Water SSB - Water Bath: Stirred (Constant) / Biological / Concentric Rings / Shaking / Refrigerated Circulating Bath - Round Bath: Clear Bottom / Immersion Series / ECO Series - Ultrasonic Bath: Digital / LCD - Dry Bath - Colony Counter - Hot Plate with Magnetic Stirrer - Anaerobic Jar: SS / Polycarbonate - Pipette & Petri Dish Sterilizing Box - Test Tube Rack - ESR Westergren & Wintrobe Apparatus - Slide Staining Stand / Rack / Trough - Slide Trays - Slide Cabinet

MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013, Tamil Nadu, India
Ph:044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz