



UV STERILISATION CABINET



EQUITRON Applications & Advantages:

- ❖ Ideal for safe keeping of sterilised items.
- ❖ Powerful ultra violet tubelight with wavelength of 253.7nm.
- ❖ Stainless steel trays of app dimension 260 x 190 x 12mm ht.
- ❖ Safety interlock on the cabinet door protects the user from accidental exposure to UV light.
- ❖ Option of wall mounting.

Model	Product	External Dim (Net) L/B/H in cm	App. Shipping Details Net / Gr. / Volume
#8105-06	UV Sterilisation Cabinet - 6 trays	34 x 30 x 37	9.0 kg / 11.0 kg / 0.086 m ³
#8105-12	UV Sterilisation Cabinet - 12 trays	34 x 30 x 67	15.5 kg / 18.0 kg / 0.136 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013. Tamil Nadu, India
Ph:044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz