



For reliable sterilisation

EQUITRON®

VERTICAL JACKETED AUTOCLAVE Standard

Applications:

- ♦ EQUITRON Standard Vertical Jacketed Autoclaves are ideal for all applications requiring, routine total destruction of all living micro-organisms.
- ♦ Ideal for Clinics, Primary Health Centres, Nursing Homes, Hospitals, Doctors and Healthcare Institutes.

EQUITRON Advantages:

- ♦ A sturdy and tension free Autoclave which can safely be operated repeatedly even in remote places.
- ♦ Comes with an additional Safety Valve for added safety.
- ♦ Water level indicator.

Construction:

- ♦ Inside chamber, jacket & outer chamber are fully stainless steel.
- ♦ Hydraulically die pressed lid from stainless steel plate.
- ♦ Flanges, fly screws, cross pins are also fully stainless steel.
- ♦ All internal joints argon welded, ground and polished to give crevice free internals.

Standard Features:

- ♦ Lid is fitted with pressure gauge, safety valve and extra safety cum manual exhaust valve.
- ♦ Moulded rubber gasket.
- ♦ Industrial grade energy efficient heater(s).

Optional Accessories:

- Stainless Steel wire mesh carrier.
- Dressing Drum.



Model	Working Chamber Size (Ø x D)	Capacity liters	OPTIONAL Carrier (Ø x D) in cm	Heater	App. Shipping Details Net / Gr. / Volume
#7427ST	30 x 50 cm (12 x 20")	35	1 of 25 x 45 - #S3084	3 KW	42 kg / 75 kg / 0.32 m ³
#7447ST	40 x 60 cm (16 x 24")	75	1 of 35 x 55 - #S3087	4 KW	85 kg / 140 kg / 0.92 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic / Jacketed, Autoclave Portable: Standard / Automatic, Autoclave SLE Series: Vacuum, Autoclave SLE series - Media Preparator - Table Top Autoclaves: Normal / Class B - UV Sterilization Cabinet - Rotary Vacuum Evaporator - Vacuum Controller - Chilled Water Circulator - PTFE Diaphragm Vacuum Pump - BOD Incubator - BSC Clean Bench - Water Still / Stone-En / Automatic Water Still - Water Bath: Stirred / Constant / Serological / Concentric Rings / Shaking / Refrigerated Circulating Bath - Round Bath: Clear Bottom / Immersion Series / ECO Series - Ultrasonic Bath: Digital / LCD - Dry Bath - Colony Counter - Hot Plate with Magnetic Stirrer - Anaerobic Jar: SS / Polycarbonate - Pipette & Petri Dish Sterilizing Box - Test Tube Rack - ESR Westergren & Wintrobe Apparatus - Slide Staining Stand / Rack / Trough - Slide Trays - Slide Cabinet

MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013, Tamil Nadu, India
Ph:044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz