



The ultimate in cleaning technology

EQUITRON[®]

ULTRASONIC BATHS

Applications:

- ❖ The **EQUITRON** Ultrasonic Baths / Cleaners deliver maximum penetration and power to address the most demanding & critical cleaning, mixing, degassing applications. Ideal for laboratory, medical, surgical, dentistry & ophthalmic instruments cleaning, printing, pharmaceutical, electroplating & allied industries.

EQUITRON Advantages:

- ❖ High frequency 53kHz for minute cavitation & deeper penetration.
- ❖ Control of Time, Temperature & Power Factor Control with heater series.
- ❖ Control of Time in the Digital Series.
- ❖ Seamless inner tanks constructed in SS 304.

Construction / Standard Features:

- ❖ Inner jointless vessel of SS 304.
- ❖ Fine tuning control of power factor between 40 & 100% with heater series.
- ❖ Inbuilt heater - range upto 60°C - with heater series.
- ❖ Built in Timer in both series.
- ❖ Non Volatile memory in both series.
- ❖ Drain Valves in models 10 litres and above.

Cleaning Solutions:

- ❖ An ideal companion for the Cleaner. Improves penetration & cleaning efficiency immensely.
- ❖ Concentrated & available in easy to use 100 ml bottles.



Model	Cap / Tank Size	Freq	U/s Power	Htr Watts	App. Shipping Details Gr.Wt. / Volume
WITH HEATER SERIES					
#8442.030.53H	03.0L / 24x14x10 cm	53 kHz	100W	75	4.0 kg / 0.009 m ³
#8442.060.53H	06.0L / 30x15x15 cm	53 kHz	180W	150	6.0 kg / 0.015 m ³
#8442.100.53H	10.0L / 30x24x15 cm	53 kHz	250W	270	8.0 kg / 0.024 m ³
#8442.150.53H	15.0L / 33x30x15 cm	53 kHz	350W	390	11.0 kg / 0.034 m ³
#8442.225.53H	22.5L / 50x30x15 cm	53 kHz	500W	520	12.0 kg / 0.051 m ³
WITHOUT HEATER SERIES					
#8422.030.53X	03.0L / 24x14x10 cm	53 kHz	100W	N.A.	4.0 kg / 0.009 m ³
#8422.060.53X	06.0L / 30x15x15 cm	53 kHz	180W	N.A.	6.0 kg / 0.015 m ³
#8422.100.53X	10.0L / 30x24x15 cm	53 kHz	250W	N.A.	8.0 kg / 0.024 m ³
#8422.150.53X	15.0L / 33x30x15 cm	53 kHz	350W	N.A.	11.0 kg / 0.034 m ³
#8422.225.53X	22.5L / 50x30x15 cm	53 kHz	500W	N.A.	12.0 kg / 0.051 m ³
SOAP CONCENTRATE FOR ULTRASONIC CLEANERS					
#S1102-05	Deconex-12, Soap Concentrate - 5 bottles of 100ml				1.5 kg / 0.005 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request. Due to continuous improvement in design, the product supplied may have modified features.



MultiLab

No.51/23, Arathoon Road, Off Appavu Lane,
Royapuram, Chennai - 600 013. Tamil Nadu, India
Ph:044-25984175, 45526026, 094441 62346
Email: multilabindia@gmail.com, www.multilab.biz



0070_0805-1-0