

ANEMOMETER DIGITAL

A1 Model AM 4201 Portable Anemometer Lutron, Taiwan make with fan type vane range 0.4 to 30m/s knots, feet and kms

A2 Model AM 4202 with temperature Lutron ,Taiwan make range 0.4 to 30m/s

A3 Model AM 4204 Hot Wire Anemometer Range 0.2 to 20 m/s and knots, feet and kms with Temp 0'C to 50'C

A4 Model AM 4206M with (Metal vane) for Velocity, Air Flow and Temperature. Lutron, Taiwan

A5 Model AM 4206 with (Plastic Vane) for Velocity, Air Flow and Temperature. Lutron, Taiwan

A6 Model AM 4210 with Temp. Lutron Make Range 0.8 to 30 m/s with RS 232 provision

A7 Model AM 4213 mini fan/vane for high temp of 80c and low velocity 0.8 to 12m/s Hot Wire

A8 Model AM 4222 Lutron, Taiwan handy with vane

A9 Cup Counter anemometer with 3 cups rotation and mechanical digit counter for wind velocity can be fixed on roof tops Deep vision Indian make as per IS 5912-1970. Gunmetal body Indian Make

A10 Digital Wind Speed Indicator system with 3 cup sensor, Range: 0.3 to 200km/h 5meter signal cable

A11 Digital Wind Speed and Wind Direction Indicator system with 3 cup sensor & wind direction sensor, Range : 0.3 to 200 km/h

A12 Model AM 4220/4221 mini vane cup Anemometer Lutron, Taiwan

A13 Anemoeter Hot Wire Equinox EQ TM-4001

A14 Model Anemometer China AVM 03 Workzone Etc Range 0-45m/Sec

A15 Mini Anemometer china Mobile Type

ANEMOMETER COMBO

A16 Model LM 8000 Taiwan 4 in1 for Velocity, Temperature, Humidity and Lux. Handy built-in vane

A17 Model LM 8010 Anemo / Temp/ Air Flow / RH / LUX 5 in 1

A18 Model LM 8102 Anemo / Sound Level / Temp / RH / LUX 5 in 1

A19 Model SP9201 / Anemo /RH / Temp / Baro / Altimete / Light/ Air Flow / Dew Point /Wet Bulb/ Wind Chill / Heat Index 11 in 1

A20 Model SP9202 / Anemo /RH / Temp. / Baro/ Altimeter/ Light /Air Flow/ Dew Point/ Wet Bulb /Wind Chill / Heat Index 11 in 1

(Pt 1000 Optional Probe Not Included)

ALTIMETER ANALOG MODEL

A21 Barigo Cat no 29 German make range 4000 meters LC 10mtrs 65mm dia (Till Stock)

A22 Range up to 4000 meters LC 20mtrs 65mm dia Japan Make

ALTIMETER DIGITAL MODEL (COMBO)

A23 Altimeter, Barometer, Thermometer and Hygrometer combine Digital Meter 4 in 1 China make

More combo Models See SL No A19 /A20

AUTOCLAVE MINI SS WITH TEFLON BEAKER

A24 Autoclave Mini Ss With Teflon Beaker 100ml

A25 Autoclave Mini Ss With Teflon Beaker 250ml

AUTOCLAVE LABORATORY

A26 Bucket type Aluminum-300x300mm (12"x12") with Pressure Gauge , Heater 20 lt

A27 Bucket type Aluminum-300x380mm (12"x15") with Pressure Gauge , Heater 20 lt

A28 Cooker Type S.Steel-300x350mm(12"x15") / Bucket Model Vertical Type with Pressure Gauge, Heater 20 lt

A29 Aluminium with S.Steel 300 x 500 Inner Aluminium Outer SS (12"x20") Double Wall Vertical Type. 40 ltrs

A30 Autoclave 250 x 450 mm fully SS with Foot Paddle lift lid operation Vertical Type W/O Accessorises 30 lt

A31 Autoclave 300 x 500 mm fully SS with Foot Paddle lift lid operation Vertical Type W/O Accessorises 40 lt

OPTIONAL ACCESSORISES FOR AUTOCLAVES SI no A23 to A 28 MODELS

A32 Digital Temp. Controller

A33 Water Level Cut Off

A34 Pressure Controller (Cut Off)

AGATE MORTAR AND PASTEL (STONE GRANITE)

A35 Agate Black Stone 2"

A36 Agate Black Stone 2.5"

A37 Agate Black Stone 3"

- A38 Agate Black Stone 3.5"
- A39 Agate Black Stone 4"
- A40 Agate Black Stone 4.5"
- A41 Agate Black Stone 5"
- A42 Agate Black Stone 5.5"
- A43 Agate Black Stone 6"
- A44 Agate Gray Stone 2"
- A45 Agate Gray Stone 2.5"
- A46 Agate Gray Stone 3"
- A47 Agate Gray Stone 3.5"
- A48 Agate Gray Stone 4"
- A49 Agate Gray Stone 4.5"
- A50 Agate Gray Stone 5"
- A51 Agate Gray Stone 5.5"
- A52 Agate Gray Stone 6"

ALUMINA WARE BOAT & CRUCIBLE

- A53 Alumina Boat - size 10x12x88
- A54 Alumina Boat - size 16x11x105
- A55 Alumina Boat - Special Size < 20mm
- A56 Alumina BoatT - special jumbo size 18X40X100 MM
- A57 Alumina Tray 24 x 63 x 122 (or) 21 x 63 x 122
- A58 Alumina Crucible 5 ml / 10 ml
- A59 Alumina Crucible 15 ml / 20 ml
- A60 Alumina Crucible 25 ml /30 ml
- A61 Alumina Crucible 50 ml
- A62 Alumina Crucible 100 ml
- A63 Alumina Crucible 150 ml
- A64 Alumina Crucible 200 ml

A65 Alumina Crucible 250 ml

A66 Alumina Crucible 500 ml

A67 Alumina Crucible 1000 ml

Allumina Lids - Prices available on request

AIR PURIFIER MUST FOR POORLY VENTILATED

ROOM / A/C ROOMS/ For SIMPLE HYGIENE A68 SANYO

Virus Washer ABC VW 24A (Till Stock Last)

ANTIBIOTIC ZONE READER

A69 Antibiotic Zone Reader Ei- Deep Vision Make Model
-1941

A70 Antibiotic ZONE Reader Duralab Digital

A71 Antibiotic ZONE Reader Duralab Analoge

AIR QUALITY ANALYZER

A72 Air quality analyzer Lutron Model AQ9901 SD , CO/Co2/
o2/ RH/ Temp. Dew Point Portable

Meter

A73 PM2.5 / PM10/Pm 1.0/ Temp /Time Air Quality
Meter YEH 310 / RT 310

A74 HCHO / PM2.5 / PM10/PM1.0/ CO2 / Temp /Time Air
Quality Meter RT 500 / YEH 500

A75 Lutron Model 1064SD /Pm2.5
/Co2/RH/Temp/Barometer

A76 CO2 / RH / Temp Model No YEH 40

A77 CO2 / RH / Temp Data logger Model YEH 40D