

CHESS CLOCKS

- C1 Chess Clock NEW Digital PS 393
- CONDUCTIVITY METER EUTECH PEN TYPE
- C2 Model ECOEC TEST LOW Range 0-1999 μ S TEMP.Display NO Resolution 10 μ S/cm Sensor Pin
- C3 Model ECOEC TEST HIGH Range 0-19.9MS TEMP.Display NO Resolution 0.1ms/cm Sensor Pin
- C4 Model EC TEST 11 Range 0-2000/ μ S & 0-20MS/cm TEMP.Display YES Resolution 10 μ S/0.1ms Sensor Pin
- C5 Model EC TEST 11 PLS Range 200,2000 μ S/cm,20MS/cm TEMP.Display YES Resolution 0.1 μ S/1 μ S/0.01ms Sensor Cup

NEW CTS EXPERT / Cond / TDS-20 ppt / Salinity / TEMP. - Refer P Multiparameter

CONDUCTIVITY METER PEN TYPE / Duralab / HM / ERMA

- C6 Model COND /TDS HM Model COM80 Range 9990 uS/5000ppm Resolution 10 μ S/10ppm TEMP.BOLD Display
- C7 Model CON Combined ERMA Model 035 Range 1990 μ S/19.99MS Resolution 1 μ S/.01ms TEMP.BOLD Display
- C8 Model CON Dist 3 HANNA Range 0-1990 μ S Resolution 1 μ S TEMP.BOLD Display
- C9 Model CON Dist 4 HANNA Range 0.01MS TO 19.99MS Resolution .01ms TEMP.BOLD Display
- C10 Model CON Tester 3 Duralab Range 1990 μ S Resolution 1 μ S Bold
- C11 Model CON Tester 4 Duralab Range 19.99MS Resolution .01ms Bold
- C12 Model CON Box CD 303 or CD 304 TESTER PEN Range 1990 μ S/19.99MS Resolution 1 μ S/.01ms Small Display
- C13 Model COND HM AP2 Range 0- 9999 μ S Resolution 1 μ S TEMP.Small Display

CONDUCTIVITY METER PORTABLE DEEP VISION /HANNA/LUTRON

- C14 HANNA HI8733 Conductivity meter Microprocessor 0 to 199.9ms (LIMITED STOCK)
- C15 LUTRON CD4302 Range 0-1.999ms, 2-19.99ms portable
- CONDUCTIVITY METER PORTABLE EUTECH
- C16 Model CON603PLUSK Parameter Conductivity Range 0 to /200MS/cm Powered by battery IP54
- C17 Model CONWP15003K Parameter Cond/TDS Range 0-200ms/200ppt Memory 150 Powered by Main/Battery

C18 Model CONWP45003k Parameter Cond/TDS/Salinity Range 200ms/200ppt/50ppt Memory 500 Powered by Main/Battery

C19 Model CONWP61043K Cond/TDS/Res/Sal Range 0 to 500MS/cm 500ppt Memory 500 Powered by Main/Battery IP 67

CONDUCTIVITY METER BENCH MODEL EUTECH

C20 Model ECCON70043S Conductivity/ATC Probe CONSEN9501D ,K=1.0 , Integral Electrode Holder & 100/240 VAC Adapter

Cond/TDS Range 0 to 200ms/cm:0 to 100 ppt Powered by Main

C21 Model ECCON270043S With 4-Cell Cond. Probe CONSEN9201D K=0.530, Integral Electrode Holder & 100/240 Adapter,

RS232 cable Cond/TDS/Salinity/Restivity and Tepm. Range 0 to 500ms/cm; 0 to 500ppt; 0 to 80 ppt,0 to 200mohms

C22 ANTECH AN 270 Printer Suitable for ECCON270043S model .

CONDUCTIVITY METER BENCH MODEL EI-DEEP VISION

C23 Model 601 Deep Vision uP Based conductivity meter 3.5 LED Display bench Model with ATC cell constant adjustment

with cell K =1.00 and temp. probe Range 0-1000 Mhos/cm in accuracy 1% with stand for cell.

C24 Model 611 Deep Vision uP Based cond. meter 3.5 LED with Manual temp.compensation & digital cell constant adj.

Range 0-1000m Mhos/cm. with cell K=1.00 w/o temp.Probe.

C25 Model 611 Advance + Deep Vision Conductivity meter with Bulit In Stirrer (Latest Model)

C26 Model 641 DV combo TDS and cond.meter table Model with both TDS and conductivity with cell K=1.00

C27 Model 1601 DV uP based 16x2 alphanumeric LCD Display showing Cond,TDS,Temp.Auto ranging,temperature co-efficient

Data Storage facility for 90 sample & external printer attachable Range Cond upto 200ms & TDS upto 200ppt in five

Selection Range with cell K=1.00 (Software For PC Communication Extra @ 6000/-)

C28 Model 181 Deep Vision combo pH and cond.and Temp meter table Model with both pH and conductivity and Temp.

C29 Model VISION - ALPHA06 Conductivity meter Economy Model (old model 611)

CONDUCTIVITY METER BENCH DURALAB

C30 Duralab Cond / TDS / Salinity / Temp. meter SCM 115 Bench Type μ P based with ATC, with Printer provision , 1000 memory

SPARE CONDUCTIVITY / TDS CELL

C31 Spare TDS/COND having cell constant at $K=1.00$ [CD-06] Banna pin [BNC against request @ Rs.175/- Extra Charges]

C32 Spare TDS/COND CELL Plastic Body $K=1.00$ Or $K = 10$ Also available

C33 Spare TDS/COND having cell constant at $K=0.50$ (appx) [CD-04]

C34 Spare TDS/COND having cell constant at $K=0.10$ (appx) [CD - 08]

C35 Spare temp probe suitable for ATC Model 101/ 601 etc.

C36 Spare metal stands to hold TDS/COND cell (or) pH Electrode REGULAR

C37 Condu.cell CONSEN 91W 2SS rings with ATC, cell constant $K=1.0$, 6 pin connector,1M cable for CON510/150/450

C38 Cond Cell CONSEN91B Suitable FOR SALT 6+ / Con 6+ / TDS 6+

C39 Cond Cell CONSEN9501D COND / ATC Electrode $K=1.0$ Suitable for Con700 Models

C40 Cond Cell CONSEN9201D 4-Cell Conductivity/ATC Electrode $K=0.530$, Graphite Sensor, Epoxy-Body (2700 series)

Calibration Buffer for conductivity Refer Page "B"

COLONY COUNTER EI-DEEPTION

(Microprocessor based Models with averaging, counting correction and 100 sample memory)

C41 Model 361 Digital Colony Counter 3 digit (999) with Magnifier & Counting Pen

C42 Model 362 Digital Colony Counter 4 digit (9999) with Magnifier & Counting Pen

C43 Model 363 Digital Colony Counter 6 digit (999999) with Magnifier & Counting Pen

C44 Model 1361 μ P Digital Colony Counter 3 digit (999) with Magnifier & Counting Pen μ P based (Micro Processor)

C45 Model 1362 μ P Digital Colony Counter 4 digit (9999) with Magnifier & Counting Pen μ P based(Micro Processor)

C46 Model 1363 μ P Digital Colony Counter 6 digit (999999) with Magnifier & Counting Pen μ P based(Micro Processor)

C47 Golden Model Colony counter

ACCESSORIES FOR COLONY COUNTER

- C48 Spare centre glass with marking
- C49 Spare Counting Pen
- C50 Spare Magnifier 100mm Dia (4") with metal frame

COLORIMETER AND HEAMOGLOBINOMETER FOR CLINICAL LAB EI DEEP VISION / ALPHA VISION COLORIMETER 8 FILTER AT 400,450,490,520,540,570,620,680 nm peak

- C51 Model 1312 Microprocessor Photo Colorimeter, 8 Glass Filters, Range 400-700nm, Auto Concentration facility % Transmission,
measurement, Absorption & K factor 116x2 line LCD Display, sample holder can accommodate Glass Test Tube,
10 mm path length cuvette, printer attachable 100 Date storage, Any Centronic
- C52 Model 1311 New μ P Photo Colorimeter 3 Digit Display with 8 glass filters AUTO ZERO
- C53 Model 312 Photo colorimeter 8 filter Digital O.D
- C54 Model 1315 Deep Vision μ P Haemoglobin meter 3 Digit LED Display with soft touch keys, auto calibration,
Direct reading in Hb units,
- C55 Model VISION ALPHA03 EI VISION Digital with 8 Filters Colorimeter (old model 312)
COUNTING CHAMBER REFER 'H' HEAMOCYTOMETER
CENTRIFUGE VISION PLUS
- C56 Alpha Vision Doctor Centrifuge 8 x 15 ml ABS body

CENTRIFUGE REMI MEDICO PLUS

- C57 Medico Plus Remi Centrifuge Max Speed 4000 RPM with Brushless Induction Motor, Frequency Drive, LED Display of
Speed and Time , Safety Lid Lock, Imbalance Detector, Dynamic Brake Supplied Without Rotors and Accessories
- C58 SM 8 - 15 MI Swing Out Rotor Suitable for 8 Tubs of 15 ml (Suppliedwith Glass tubes.
CENTRIFUGE MINI / MICRO REMI Revised On 15.06.2021
- C59 RM 02 PLUS Remi Mini Plus Centrifuge 6000 RPM With Speed Indicator Capacity 8 x 1.5 ml
- C60 RM 03 PLUS Remi Micro Centrifuge 15000 RPM LCD Display with 12x1.5ml / 2 ml and PCR strip Rotor (2 x 8 x 0.2) ml

CENTRIFUGE REMI

C61 C-852 Medical Centrifuge Buckets 4 Capacity 15ml RPM 3500 NO Display

C62 C-854/4 Clinical Centrifuge Bucket 4 Capacity 15ml RPM 3500 NO Display

C63 C-854/6 Clinical Centrifuge Bucket 6 Capacity 15ml RPM 3500 NO Display

C64 C-854/8 Clinical Centrifuge Bucket 8 Capacity 15ml RPM 3500 NO Display

C65 R-303 Doctors Centrifuge Bucket 8 Capacity 15ml RPM 3800 NO Display

C66 R-304 Doctors Centrifuge Bucket 24 Capacity 5ml RPM 7000 NO Display

C67 R-4C W/Timer Centrifuge Without Head RPM 4400 WITH Display

C68 Head R-41 Spare Bucket 8 (SWING) Capacity 15ml RPM 3500, Suitable for R-4C

C69 Head R-42 Spare Bucket 4 (SWING) Capacity 50ml RPM 3500 Suitable for R-4C

C70 Head R-43 Spare Bucket 8 (ANGLE) Capacity 15ml RPM 4400 Suitable for R-4C

C71 Head R-44 Spare Bucket 4 (ANGLE) Capacity 50ml RPM 4200 Suitable for R-4C

C72 R-8C Centrifuge Without Head RPM 5250 WITH Display

C73 R8C BL Centrifuge Without Head RPM 6000 WITH Display

C74 Head R-81 Spare Bucket 16 (swing) Capacity 15ml RPM 4000 Suitable for R-8C

C75 Head R-81A Spare Bucket 12(swing) Capacity 15ml RPM 4000 Suitable for R-8C

C76 Head R-81B Spare Bucket 8(swing) Capacity 15ml RPM 4200 Suitable for R-8C

C77 Head R-82A Spare Bucket 4 (swing) Capacity 50ml RPM 4000 Suitable for R-8C

C78 Head R-83A Spare Bucket 12 (Angle) Capacity 15ml RPM 5250 Suitable for R-8C

C79 Head R-84A Spare Bucket 4 (Angle) Capacity 50ml RPM 5250 Suitable for R-8C

C80 Head R-88 Spare Angle Head 4 x 100 ml with Polypropylene tube For R-8C BL

C81 Head R-89 Spare Angle Head 24 x 15 ml with Polypropylene tube For R-8C BL

C82 Head R-90 Spare 4 x 100 ml swing out with glass tube & Biosafety cap for R-8CBL

C83 Head R-82 Spare Bucket 6 (swing) Capacity 50ml Suitable for R-8C

C84 Head R-83 Spare Bucket 16 (Angle) Capacity 15ml Suitable for R-8C

C85 RM-12C W/timer Micro Centrifuge w/o Head RPM 16000 with Display

C86 Head RM-1210 Spare Bucket 24 Capillary Capacity 75mmx1mm RPM 12000 Suitable for RM - 12C

C87 Head RM-1214 Spare Bucket 10 (Angle) Capacity 2.0ml RPM 16000 Suitable for RM - 12C

C88 Head RM-1215 Spare Bucket 24 (Angle) Capacity 1.5ml RPM 11850

- C89 Head RM-1213 Spare
- C90 Head RM 1216 Spare Bucket 8 Capacity 5ml
- C91 R-8K Centrifuge For Milk & Oil W /O Rotor with Display
- C92 Head R-86 Spare Rotor for above suitable for R-8k Capacity 4 x 50 ml Head
- C93 R-8M PLUS Centrifuge (Max 6000 RPM) With Brushless Induction Motor without Head
- C94 R-8C PLUS Centrifuge (Max 6000 RPM) With Brushless Induction Motor without Head
- C95 Head S 24 -15 Swing Out Rotor 24 x 15ml (Supplied with Glass Tubes)
- C96 Head S 8 -15 Swing Out Rotor 8 x15ml (Supplied with Glass Tubes)
- C97 Head AR24 - 15 Angle Rotor 24 x 15 ml (Supplied with Plastic Tubes)
- C98 Head AR6 - 50 Angle Rotor 6 x 50 ml (Supplied with Plastic Tubes)

CENTRIGUGE REMI NEYA SERIES

- C99 NEYA 12 Bench Top Centrifuge, (max. speed 16000 rpm) LCD display of speed and time, imbalance detector,

dynamic brake safety lid lock, 10 acceleration / deceleration profiles, 10 programs and rotor identification,

supplied without rotors and accessories
- C100 NEYA 16 Refrigerated Bench Top Centrifuge, (max. speed 16000 rpm) LCD display of speed and time,

imbalance detector, dynamic brake, safety lid lock, 10 acceleration / deceleration profiles,

10 programs and rotor identification, supplied without rotors and accessories
- C101 Head A 24-2 /24 x 2ml Angle Rotor for 2ml Micro Tubes with Reducers RE 2-05 (for 0.5ml) & RE 2-02 (for 0.2ml)
- C102 Head A 24-15/ 24x15ml Angle Rotor (supplied with plastic tubes)
- C103 S 4-175 4 place Swing out Rotor
- C104 S 6-96 MP 2 place Swing out Rotor for 6 micro plates of 96 wells

CYCLOMIXER / VORTEX MIXER REMI AND REMI PLUS SERIES

- C105 REMI Model CM101 with single and Multi tube Head
- C106 CM-101 Plus Vortex Mixer with brushless motor & suitable for Interchangeable mixing Heads

CENTRIFUGE HIGH SPEED & COOLING & REFRIGERATED NORMAL& PLUS SERIES REMI

- C107 R-24 Revolutionary High Speed Centrifuge Without Rotor Head And Accessories
- C108 R-242 ANGLE Head Rotor With PP Tubes 6X50ML Suitable For R-24 And C-24
- C109 CM-12 Plus Refrigerated Compact Bench Top Cooling Micro Centrifuges w/o rotor Heads and accs -20'c up to 5 ml)
- C110 RM-1214M 10 x 2 ml Angle Head (with Polycarbonate Tubes)
- C111 RM-1215M 24 x 1.5 ml Angle Head (with Polycarbonate Tubes)
- C112 RM-1216M 8 x 5 ml Angle Head (with Polycarbonate Tubes)
- C113 RM-1217M Angle Head suitable for PCR Strip Tubes - Capacity 4 strips of 8 x 0.2 ml
- C114 VS-02 Servo controlled Voltage Stabilizer
- C115 C-23 Plus HIGH Volume Low Speed Table Top Cooling Centrifuge Max 7000 RPM / 7560 RCF, LED Display w/o rotor Heads
- C116 CPR-23 Plus HIGH Volume Low Speed Table Top Cooling Max 7000 RPM / 7400 RCF, LCD Display W/O rotor Heads
- C117 R-238M Complete with Polypropylene Tubes (Capacity 12 x 50 ML)
- C118 R-239M Complete with Polypropylene Tubes (Capacity 36 x 15 ML)
- C119 VS-03 Servo controlled Voltage Stabilizer
- C120 C-24 Plus Low Volume High Speed Cooling MAX 20000 RPM / 37570 RCF, LED Display W/O rotor Heads
- C121 CPR-24 Plus Low Volume High Speed Cooling, MAX 20000 RPM / 37570 RCF, LCD Display W/O rotor Heads
- C122 C-30 Plus HIGH Volume High Speed FLOOR Model COOLING MAX 20000 RPM / 37570 RCF, LED Display, W/O rotor
- C123 CPR-30 Plus HIGH Volume High Speed FLOOR Model COOLING MAX 20000 RPM / 37570 RCF LCD Display w/o Rotor

LCD Display W/O rotor Heads

- C124 R-243M Complete with Polypropylene Tubes (Capacity 8 x 25 ML)
- C125 R-244M Complete with Polypropylene Tubes (Capacity 12 x 15 ML)
- C126 R-247M Complete with Polypropylene Tubes (Capacity 8 x 50 ML)
- C127 R-248M Complete with Polypropylene Tubes (Capacity 24 x 1.5 ML)
- C128 VS-03 Servo controlled Voltage Stabilizer

COLD CABINETS (LABORATORY REFRIGERATOR) REMI

- C129 CC-6 Plus Capacity 200 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02)
LED
- C130 CC-6 Plus Capacity 200 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02)
LCD
- C131 CC-10 Plus Capacity 280 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02)
LED
- C132 CC-10 Plus Capacity 280 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02)
LCD
- C133 CC-12 Plus Capacity 370 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02)
LED
- C134 CC-12 Plus Capacity 370 liters, Temperature Range 2°C to 12°C Microprocessor based (VS-02)
LCD
- C135 VS-02 Servo controlled Voltage Stabilizer

The high end REMI Items will be Supplied Against Order With Advance . Delivery 3 - 4 Weeks

CARBON MONOXIDE METER (CO)

- C136 Carbon Monoxide Range 0 - 1000 ppm with Temp - 0 - 50°C Model Lutron GCO 2008
CO2 METER
- C137 CO2 Meter Model GCH-2018 For Measuring Carbon Dioxide Range 0-3000ppm With Alarms And
Warning. RH 10 to 95%,
temp, 0 to 50c and dew point Lutron Taiwan Make
- C138 CO2 meter GC- 2028 With alarm Setting Lutron Taiwan Make
- C139 CO2 / RH / Temp Model No YEH 40
- C140 CO2 / RH / Temp Data logger Model YEH 40D
CO2 AIR QUALITY METER Refer Page "A"
CAPILLARY TUBES TOP TECH
HAEMOTOCRIT CAPILLARY
- C141 75mm Non Heparinised Capillary 1 / 100Nos
- C142 75mm Non Heparinised Capillary 1 / 250Nos
- C143 75mm Heparinised Capillary 1 / 100Nos
- C144 75mm Heparinised Capillary 1 / 250Nos

C145 125mm Non Heparinised Capillary BGA 1 / 100Nos

C146 125mm Heparinised Capillary BGA 1 / 100Nos

MELTING POINT CAPILLARY

C147 Both End Open 3" Long 1 / 100Nos

C148 Both End Open 3" Long 1 / 250Nos

C149 Both End Open 4" Long 1 / 100Nos

C150 Both End Open 4" Long 1 / 250Nos

C151 One End Sealed 3" Long 1 / 100Nos

C152 One End Sealed 3" Long 1 / 250Nos

C153 One End Sealed 4" Long 1 / 100Nos

C154 One End Sealed 4" Long 1 / 250Nos

CAPILLARY GRADUATED PER PACK CAPILLARIES (Pkg of 250 nos) TOP TECH

C155 5ul 5 mark with colour code

C156 10ul 2 mark with colour code

C157 5/ 10/ 20/ 25/ 50/ 100ul Graduated Capillary 1 mark with colour code

C158 5/ 10/ 20/ 50/ 100ul Graduated Capillary without colour code 90mm

C159 TLC Capillary thin without graduation TOPTECH (1 x 250)

COMPARATOR WITHOUT COLOUR DISCS LOVIBOND

C160 COMPARATOR 2000+ with 1 pair of square 13.5 mm cells LOVIBOND Code No 142000 Or PALINTEST

C161 Color Disc Lovibond PH / SILICA / CHLORINE etc. Each One

C162 Spare 13.5mm, 10ml moulded cells for comparator with stoppers PLASTIC (One Pair) LOVIBOND

C163 DPD Tablets No.1/2/3 or 4 LOVIBOND Pkt of 250 Tablet

C164 Nesselizer 2150 with stand (without Day Light) Code No 172150

COMPARATOR & DISCS INDIAN

C165 COMPARATOR with 2 Color disc for Lub Oil ASTM D 1500

C166 NESSELIZER INDIAN WITH 2 COLOR DISC FOR LUB OIL ASTM D 1500

C167 COMPARATOR with 1 pair of square 13.5 mm cells "INDIAN"

C168 NESSELIZER INDIAN WITH SAMPLE TUBES

- C169 Spare colour disc price code A CHLORINE / pH / UNIVERSAL INDIAN each Disk
- C170 Spare colour disc price code B - HAZEN INDIAN
- C171 Spare colour disc price code C INDIAN
- C172 Spare 13.5mm, 10ml moulded cells for comparator with stoppers (One Pair) INDIAN
- C173 DPD Tables No.1/2/3 or 4 (Price per Each Tablet) LOVIBOND

CHECKIT COMPARATOR TEST KIT PALINTEST (Comparator Disc and Reagents) Plastic Color Discs

- C174 CKH 1001 Contour Comparator kit, Free and Total Chlorine Range 0-5mg/l Cl₂
- C175 CKH 1192 Contour Comparator kit, Alkavis Alkalinity Range 0-250mg/l CaCO₃
- C176 CKH 1152 Contour Comparator kit, Ammonia Range 0-1mg/l N
- C177 CKH 1179 Contour Comparator kit, Fluoride Range 0-1.5mg/l F
- C178 CKH 1114 Contour Comparator kit, Phosphate HR Range 0-100mg/l PO₄
- C179 CKH 1181 Contour Comparator kit, Silica Range 0-4mg/l SiO₂
- C180 CKH 1292 Contour Comparator Kit, Iron Range 0.10mg/l Fe

Chemical Test Kit Refer " Page T "

COMBUSTION TUBE

- C181 17x22x600 kumar Make
- C182 19x24x600 Kumar Make

COMBUSTION BOAT

- C183 KUMAR size 10x12x88mm (pkt of 150nos)- - Price per piece
- C184 Ordinary quality size 10x12x88mm (pkt of 100nos) Price per piece

HAND CENTRIFUGE

- C185 2 TUBE
- C186 4 TUBE

DP PENCIL / LOOP HOLDER / TEST TUBE HOLDER

- C187 Test Tube Holder
- C188 DP Pencil
- C189 Loop Holder

COD MEASUREMENT SYSTEM

C190 COD Measurement system COD set up MD 200 COD VARIO 10 with reagents consisting of COD photometer, LOVIBOND

COD reactor RD 125Code 2418940, COD tube 0 - 150 mg/l and 0-1500 mg/l cuvette rack, manual code 2892602

C191 COD ECO8 F101A0127 VELP Make Digester with 25 Vario Tube Reagent 0 - 2000 mg/l with Multiparameter Photometer

PALINTEST PTH7100

C192 Lovibond COD Vario tubes 0-150 mg/l or 0-1500 mg/l /15000 mg/l (Reagents) 2420710 / 2420711/2420712 each pack of 25

C193 COD Digester 16 test tube VELP Make ECO-16

C194 COD Digester kit single tube small size Digester with reagents and vials Indian

COD DIGESTER INDIAN DEEPCISION

C195 COD Digester with digital temp controller & timer along with glass parts Indian make 15 Samples

C196 COD Digester with digital temp controller & timer along with glass parts Indian make 6 Samples

CHLORINE METER DIGITAL (Also See ION Meters in I for other Parameters)

C197 Chlorine Meter (Free & Total _ Lutron CL 2006 Range 0 - 3.5 mg/l with 10 test reagent free and total STD included

C198 Chlorine Meter Portable Digital Aquasol Model AP-CL-01 Range 0 - 5 ppm, Resolution 0.01 ppm

C199 Chlorine Handy Photometer for Free Chlorine range 0 to 2.5ppm Milwaukee MW10

C200 Spare Reagent (250 Test) Free Chlorine Cat No AP-CL-FC1

C201 Spare Reagent (250 Test) Total Chlorine Cat No AP-CL-TC1

COATING THICKNESS GAUGE DIGITAL

C202 Model M1Coating thickness gauge AA make digital Range 2000 microns Ferrous to findout nonmagnetic Coating

on a magnetic surface

C203 Model C1Coating thickness gauge AA make digital Range 2000 microns Non Ferrous to findout nonmagnetic Coating

on a Nonmagnetic surface

C204 Model M&C 1 AA make digital Range 2000 microns to find out non magnetic coating on a magnetic surface and conductive

coating on a not magnetic surface

C205 Coating thickness gauge with inbuilt sensors F/NF with data output and software, China make DURALAB / AccuPlus

C206 Workzone DFT 222 Coating thickness gauge F/NF with Inbuilt Sensors with data output and software, China make

C207 Metravi CTG-01 Coating thickness gauge with inbuilt sensors for ferrous and non ferrous metals with data output & software

C208 Workzone CoatPro Coating thickness gauge for Ferrous and Non Ferrous metals wire probe

C209 Workzone DFT 456 Coating thickness gauge with separate SS probe for both Ferrous and Non Ferrous

C210 Workzone COATFNF Coating thickness gauge with Sensor with bluetooth China make Range 0-1500microns

COATING THICKNESS GAUGE SPARES

C211 Spare standard coating foils

C212 Spare standard base for zero setting

COATING THICKNESS GAUGE MAGNETIC BASED MANUAL MODELS

C213 Coating thickness Gauge 0-1000 microns Elcometer for ferrous materials

C214 Coating Thickness Pencil Type manual upto 600 microns

C15 Coating Thickness Pencil with Cursor upto 200 microns

Carbon Residue Apparatus - " Please Refer " P " Petroloum Instruments.

CHRONO METER

C216 Chronometer Clock Barigo German Cat No 1816

CLAMP METER

C217 Clamp Meter Work Zone Model 3600 Ac/Dc Max 1000 V Dc and 750 V AC