

DISSOLVED OXYGEN METER PEN TYPE LUTRON

- D1 Model PDO 519 Lutron DO meter Pen type 0-20 ppm, Taiwan with built in probe
- D2 Model PDO 520 Lutron DO meter Pen type 0-20 ppm, Taiwan with Detachable probe

D.O METER PORTABLE MODEL - LUTRON

- D3 Model 5509 Lutron DO Meter, 0 to 20mg/L manual calibration
- D4 Model 5510 Lutron DO Meter 0 to 20mg/L measure oxygen in air / auto calibration
- D5 Model 5510 HA Lutron DO meter 0 - 20 mg/l / Oxygen in Air - 0 - 100 %

DISSOLVED OXYGEN METER PORTABLE MODEL- EUTECH and LOVIBOND

- D6 Model DO602K DO Range 0 to 19.9 mg/l; 0 TO 199.9% IP Rating IP 54
- D7 Model SD310 Oxi Cat No. 724650 DO meter Portable Lovibond with probe 1.5m cable (Only One Unit Left)

D.O METER BENCH MODEL- EUTECH

- D8 Model ECDO70042S DO/Temp RANGE 0 to 30 mg/l NO Memory
- D9 Model ECDO270042 DO/BOD/Temp. RANGE 0 to 50 mg/l Memory 500 Output RS232
- D10 HANNA DO / BOD Meter Model HI 5421-02 (ONLY ONE UNIT LEFT)

D.O METER - ACCESSORIES

- D11 Spare DO probe for Lutron 5509/5510
- D12 Spare DO membrane for Lutron DO (Pack of 2 Replaceable membrane)
- D13 Spare Electrolyte for Lutron DO meters 30 ml bottle

DIFFRACTION GRATING

- D14 Diffraction grating 15000 Lines/ inch OR 7500 lines OR 2500 lines Imported China Make
- D15 Diffraction grating 15000 lines/inch Indian
- D16 Diffraction grating 2500 Lines/inch Indian for laser
- D17 Diffraction grating 7500 Lines/inch Indian
- D18 Diffraction 3 in One 15000 / 7500 / 2500 lines / inch Indian
DEAN AND STARK APPARATUS (Also Refer 'P' Petroleum Testing Apparatus)
- D19 To know the moisturisation in Lubricating oils as per IP 74 ASTM D 95 IS 1448

DEHUMIDIFIER

- D20 DEHUMIDIFIER New Model ADVANCE A3DF200 (Till Stock Last)

DRY BATH (INCUBATOR)

D21 Dry Bath Incubator Deep Vision Model 352 ambient +5'C to 100'C with fixed tube trays size 12 mm x 24 holes

DEEP FREEZERS(VERTICAL) REMI PLUS SERIES (Against Order with Advance only) Revised On 15.06.2021

D22 RQV-200 Plus Capacity 200 liters, Temperature up to -20°C Microprocessor (VS 02) LED

D23 RQV-200 Plus Capacity 200 liters, Temperature up to -20°C Microprocessor (VS 02) LCD

D24 RQV-300 Plus Capacity 300 liters, Temperature up to -20°C Microprocessor (VS 02) LED

D25 RQV-300 Plus Capacity 300 liters, Temperature up to -20°C Microprocessor (VS 02) LCD

D26 RQV-400 Plus Capacity 400 liters, Temperature up to -20°C Microprocessor (VS 03) LED

D27 RQV-400 Plus Capacity 400 liters, Temperature up to -20°C Microprocessor (VS 03) LCD

D28 VS-03 Servo controlled Voltage Stabilizer

D29 VS-02 Servo controlled Voltage Stabilizer

ULTRA DEEP FREEZERS (VERTICAL) REMI SERIES (Against Order with Advance)

D30 ULT-90 Capacity 90 liters Temp low upto -86c microprocessor control lcd version(VS 05)

D31 ULT-185 Capacity 185 liters Temp low upto -86c microprocessor control lcd version(VS 06)

D32 ULT-360 Capacity 360 liters Temp low upto -86c microprocessor control lcd version(VS 07)

D33 VS-05 Servo controlled Voltage Stabilizer

D34 VS-06 Servo controlled Voltage Stabilizer

D35 VS-07 Servo controlled Voltage Stabilizer

DISTANCE METER - (LAZER TYPE)

D36 Work Zone / Equal China Make LDM - 60 (60 Meter)

D37 Work Zone / Equal China Make LDM - 80 (80 Meter)

D38 Work Zone / Equal China Make LDM - 100 (100 Meter)