

HOT AIR OVEN DEEP VISION

Double walled Outer body made of heavy gauge mild steel. Sufficient glass wool insulation. Complete with pilot

lamps, rotary switch to work on 230Volts AC , Provided with expanded metal Shelves, plug, ord and connector.

Thermostat control temp Ambient Range to 250'C Inner SS With Outer MS Power Coated Body

H1 Chamber size 300x300x300mm 12"x12"x12" capacity 28 Ltr,750 Rating 1 Shelves Analog Regulator Temp.

H2 12"x12"x12" with Digital Temperature Controller

H3 12"x12"x12" with Digital Temperature Controller and Fan for Circulating air

H4 Chamber Size 355x355x355mm 14"x14"x14" capacity 45 Ltr,1000 Rating 2 Shelves Analog Regulator Temp.

H5 14"x14"x14" with Digital Temperature Controller

H6 14"x14"x14" with Digital Temperature Controller and Fan for Circulating air

H7 Chamber size 455x455x455mm 18"x18"x18" capacity 95 Ltr, 1500 Rating, 2 Shelves, Analog Regulator Temp.

H8 18"x18"x18 with Digital Temperature Controller

H9 18"x18"x18 with Digital Temperature Controller and Fan for Circulating air

H10 Inner Chamber 605x605x605mm 24"x24"x24" capacity 224 Ltr, 2200 Rating, 3 Shelves Analog Regulator Temp.

H11 24"x24"x24" with Digital Temperature Controller

H12 24"x24"x24" with Digital Temperature Controller and Fan for Circulating air

H13 Optional Accessories PID Digital Temp Controller for Better Accuracy

HOT AIR OVEN

H14 RDHO-50 DryHot Air Oven (50 liters) with temperature setting from RT + 10°C up to +300°C

H15 RDHO-80 DryHot Air Oven (80 Liters) with temperature setting from RT + 10°C up to +300°C

HOT PLATE DEEP VISION

H16 Hot Plate8" Dia (1000W) M.S. (Energy Regulator)

H17 Hot Plate9" Dia (1000W) M.S. (Energy Regulator)

H18 Hot Plate12" Dia (1500W) M.S. (Energy Regulator)

H19 Hot Plate10x12 (1000W) M.S. (Energy Regulator)

- H20 Hot Plate10x16 (1500W) M.S. (Energy Regulator)
- H21 Hot Plate12x18 (2000W) M.S. (Energy Regulator)
- H22 Hot Plate18x24 (4000W) M.S. (Energy Regulator)
- H23 Hot Plate 8" Digital Controller MS BODY
- H24 Hot Plate12X18" Digital Controller M.S Digital
- H25 Hot Plate7" [ORDINARY] M.S Body Commercial
- H26 Hot Plate 8" Dia (1000W) S.S. (Energy Regulator)
- H27 Hot Plate9" Dia (1000W) S.S. (Energy Regulator)
- H28 Hot Plate12" Dia (1500W) S.S. (Energy Regulator)
- H29 Hot Plate10x12 (1000W) S.S. (Energy Regulator)
- H30 Hot Plate10x16 (1500W) S.S. (Energy Regulator)
- H31 Hot Plate12x18 (2000W) S.S. (Energy Regulator)
- H32 Hot Plate12x24 (4000W) S.S. (Energy Regulator)
- H33 Hot Plate 8" Digital Controller SS body
- H34 Hot Plate12x18" SS Top Digital Controller
- H35 Digital controller for any one of the above extra charge

HEATING MANTELS DEEP VISION

- H36 Heating Mantels Capacity 50ml
- H37 Heating Mantels Capacity 100ml
- H38 Heating Mantels Capacity 250ml
- H39 Heating Mantels Capacity 500ml
- H40 Heating Mantels Capacity 1000ml
- H41 Heating Mantels Capacity 2000ml
- H42 Heating Mantels Capacity 3000ml
- H43 Heating Mantels Capacity 5000ml
- H44 Heating Mantels Capacity 10000ml
- H45 Heating Mantels Capacity 20000ml
- H46 Spare Mantel 50ml
- H47 Spare Mantel 100ml

- H48 Spare Mantel 250ml
- H49 Spare Mantel 500ml
- H50 Spare Mantel 1000ml
- H51 Spare Mantel 2000ml
- H52 Spare Mantel 3000ml
- H53 Spare Mantel 5000ml
- H54 Spare Mantel 10000ml
- H55 Spare Mantel 20000ml

HEATING MANTEL WITH MAGNETIC STIRRER REMI

- H56 Rota Mantel REMI 1RML Capacity 1000ml
- H57 Rota Mantel REMI 2RML Capacity 2000ml

HEGMAN GAUGE as per IS 101 made of Stainless Steel

- H58 0-100microns / 0-50microns / 0-25microns / 0-10 Microns- Any One
- H59 Hegman Gauge Specified Range With Certificate

HARDNESS TESTER FOR RUBBER

- H60 Rubber hardness tester Shore A KORI, Japan make Dial type
- H61 Rubber Hardness Analog Shore A Traingular STI Make
- H62 Rubber Hardness Analog Shore D Traingular STI Make
- H63 Rubber Hardness Tester ISH-SAM Dial Shore A Model INSIZE
- H64 Rubber Hardness Tester ISH-SDM Dial Shore D Model INSIZE
- H65 Rubber hardness Tester Shore A Analog China
- H66 Rubber Hardness Tester Shore A Digital China/ or / Shore D digital China make

HYDROMETER DENSITY LIEMCO ISI - AMBALA MAKE

Hydrometer Density for petroleum calibrated at 15c I.S.I. Mark 3104

- H67 DN Range 0.600 to 0.650 or 1.000 to 1.050 LIEMCO, AMB M-50(0.001)
 - H68 DN Range 0.600 to 0.650 or 1.000 to 1.050 LEIMCO/LIEMCO L-50(0.0005)
 - H69 DN Range 0.650 to 0.700 Or 0.750 to 0.800 LIEMCO, AMB M-50(0.001)
 - H70 DN Range 0.650 to 0.700 Or 0.750 to 0.800 LIEMCO,AMB L-50(0.0005)
 - H71 DN Range 0.700 to 0.750 for PETROL LIEMCO, AMB M-50(0.001)
 - H72 DN Range 0.700 to 0.750 for PETROL LIEMCO,AMB L-50(0.0005)
 - H73 DN Range 0.800 to 0.850 for Diesel LIEMCO, AMB M-50(0.001)
 - H74 DN Range 0.800 to 0.850 for Diesel LIEMCO,AMB L-50(0.0005)
 - H75 DN Range 0.850 to 0.900 LIEMCO, AMB M-50(0.001)
 - H76 DN Range 0.850 to 0.900 LIEMCO,AMB L-50(0.0005)
 - H77 DN Range 0.900 to 0.950 or 0.950 to 1.000 LIEMCO, AMB M-50(0.001)
 - H78 DN Range 0.900 to 0.950 or 0.950 to 1.000 LIEMCO,AMB L-50(0.0005)
 - H79 Calibration Certificate by private for Hydrometer with Traceability Xerox
 - H80 Glass Cylinder 150ml with stopper A class Omson make
 - H81 Glass Cylinder 100ml with stopper A class Omson make
 - H82 Glass Graduated Pipette 10ml with A class Omson make
 - H83 Glass Cylinder with / without markings for Hydrometers borosilicate clear glass 500 ml
 - H84 Water Paste OR Oil Paste for Density Hydrometer OLEODIS (Mumbai)
 - H85 Water Paste OR Oil Paste for Cheap [ambala]
 - H86 Glass Measuring Jar for Petrol Bunk 5 LTR (5000 ml) certified
 - H87 For Petrol bunk Hydrometer kit (Consist of 2 Hydrometer / Water Paste / Oil Paste / Measuring Jar / Filter Paper)
 - H88 Sample Bottole Continer
 - H89 Brass Measuring Jar for Petrol Pump use 5 LTR
- PETROL DENSITY HYDROMETER CERTIFIED
- H90 Liemco M.50 Range DN 0.700 to 0.750 certified (Div 0.001)
 - H91 Liemco L.50 Range DN 0.700 to 0.750 certified (Div 0.0005)

DIESEL DENSITY HYDROMETER CERTIFIED

- H92 Liemco M.50 Range DN 0.800 to 0.850 certified (Div 0.001)
- H93 Liemco L.50 Range DN 0.800 to 0.850 certified (Div 0.0005)
- HYDROMETERS (IMPORTED)
- H94 Alcohol Meter 0 - 100% GAY LUCAS GERMAN
- H95 Specific Gravity Hydrometer Difference Range GERMAN 700-1000/1000-2000/others Any one
- H96 Hydro With Thermometer brix Range GERMAN - Range 0-10/10-20/20-30 (or) 0-30 Any one (Delivery by hand only from our Office)
- H97 Brix Hydrometer 0-10/10-20/20-30 GERMAN - Any One
- H98 Baume Hydrometer 0-10/10-20/20-30 & 0-30 GERMAN Any One
- H99 Sikes Hydrometer Zeal England Range 0-20A/0-20/20-40/40-60/60-80/80-100 Any One [Hand Delivery] Any One
- H100 Zeal Hydrometer 1000-2000/700-1000/1100-1300/ Alcohol Hydrometer Gay Lucas 0-100% Any one

HYDROMETER LIEMCO, AMBALA

Specific Gravity, Baume, Brix, Soil, Salino, Alcohol, Sikes, Twaddel Indian (Leimco / Liemco Brands) Any Available One

- H101 Hydrometer Range SP,Gravity 1.000-2.000 (0-70) LEIMCO
- H102 Hydrometer Range SP,Gravity .700 - 1.000 LEIMCO
- H103 Hydrometer Range SP,Gravity 1.000 - 1.500 LEIMCO
- H104 Hydrometer Range SP,Gravity 1.500 - 2.000 LEIMCO
- H105 Hydrometer 0.200 Range 0.600-.800 LEIMCO
- H106 Hydrometer 0.200 Range 0.800 - 1.000 LEIMCO
- H107 Hydrometer 0.200 Range 1.000 - 1.200 LEIMCO or 0.800 - 1.000
- H108 Hydrometer 0.200 Range SP,Gravity 1.200 - 1.400 LEIMCO
- H109 Hydrometer 0.200 Range SP,Gravity 1.400 - 1.600 LEIMCO
- H110 Hydrometer 0.200 Range SP,Gravity 1.600 - 1.800 LEIMCO
- H111 Hydrometer 0.200 Range SP,Gravity 1.800 - 2.000 LEIMCO
- H112 Hydrometer 0.100 Range 0.600 -.700
- H113 Hydrometer 0.100 Range 0.700 to 0.800/ 0.001

- H114 Hydrometer 0.100 Range 0.800 - 0.900 / 0.900 - 1.000
- H115 Hydrometer 0.100 Range 1.000 - 1.100
- H116 Hydrometer 0.100 Range 1.100 - 1.200
- H117 Hydrometer 0.100 Range 1.200 - 1.300
- H118 Hydrometer 0.100 Range 1.300 - 1.400
- H119 Hydrometer 0.100 Range 1.400-1.500
- H120 Hydrometer 0.100 Range 1.500-1.600
- H121 Hydrometer 0.100 Range 1.600 - 1.700
- H122 Hydrometer 0.100 Range 1.700 - 1.800
- H123 Hydrometer 0.100 Range 1.800 - 1.900
- H124 Hydrometer 0.100 Range 1.900 - 2.000
- H125 Hydrometer 0.050 Range 0.600 - 0.650
- H126 Hydrometer 0.050 Range 0.650 - 0.700
- H127 Hydrometer 0.050 Range 0.700 - 0.750
- H128 Hydrometer 0.050 Range 0.750 - 0.800
- H129 Hydrometer 0.050 Range 0.800 - 0.850
- H130 Hydrometer 0.050 Range 0.850 - 0.900
- H131 Hydrometer 0.050 Range 0.900 - 0.950
- H132 Hydrometer 0.050 Range 0.950 - 1.000
- H133 Hydrometer 0.050 Range 1.000 - 1.050
- H134 Hydrometer 0.050 Range 1.050 - 1.100
- H135 Hydrometer 0.050 Range 1.100 - 1.150
- H136 Hydrometer 0.050 Range 1.150 - 1.200
- H137 Hydrometer 0.050 Range 1.200 - 1.250
- H138 Hydrometer 0.050 Range 1.250 - 1.300

SALINO Hydrometer

- H139 SALINO Hydrometer 0-100%TELTRU

BAUME Hydrometer

- H140 Hydrometer 0 - 10
- H141 Hydrometer 10 - 20
- H142 Hydrometer 20-30
- H143 Hydrometer 30-40
- H144 Hydrometer 40-50
- H145 Hydrometer 50-60
- H146 Hydrometer 60-70
- H147 Hydrometer 0-20
- H148 Hydrometer 20-40
- H149 Hydrometer 40-60
- H150 Hydrometer 0-30
- H151 Hydrometer 30-60
- H152 Hydrometer 0-40
- H153 Hydrometer 0-50
- H154 Hydrometer 0-70

BRIX Hydrometer

- H155 BRIX Hydrometer 0-10
- H156 BRIX Hydrometer 10-20
- H157 BRIX Hydrometer 20-30
- H158 BRIX Hydrometer 0-30
- H159 BRIX Hydrometer 30-40
- H160 BRIX Hydrometer 30 - 60
- H161 BRIX Hydrometer 40 - 50
- H162 BRIX Hydrometer 50 - 60
- H163 BRIX Hydrometer 60 - 70
- H164 BRIX Hydrometer 60 - 90
- H165 BRIX Hydrometer 90-100

BRIX With Thermo

- H166 BRIX With Thermo 0 - 10
- H167 BRIX With Thermo 10 - 20
- H168 BRIX With Thermo 20 - 30
- H169 BRIX With Thermo 0 - 30
- H170 BRIX With Thermo 30 - 60
- H171 BRIX With Thermo 60 - 90

SOIL Hydrometer

- H172 SOIL Hydrometer 0.995 TO 1.030 LEIMCO
- H173 TELTRU

Alcohol hydro- Meter

- H174 Alcohol Meter 0 - 100% (Gay Lucas)
- H175 Alcohol Meter 90 - 100%
- H176 Alcohol Meter 94 - 104
- H177 Alcohol Meter 0-100%TELTRU

Sikes

- H178 Sikes 0 -20A
- H179 Sikes 0 - 20
- H180 Sikes 20 - 40
- H181 Sikes 40 - 60
- H182 Sikes 60 - 80
- H183 Sikes 80 - 100
- H184 Sikes TABLE BOOK

Barko Meter

- H185 Barko Meter 0 - 120

Twaddle

- H186 Twaddle No.1 (0 - 26)
- H187 Twaddle No.2 (24 - 50)
- H188 Twaddle No.3 (48 - 74)

- H189 Twaddle No.4 (72 - 104)
- H190 Twaddle No.5 (102 - 138)
- H191 Twaddle No.6 (136 - 174)
- H192 Twaddle No 7 (160-200)

Battery Hydrometer

- H193 Battery Hydrometer Big Leimco-Range 1.100 to 1.300 specific Gravity
- H194 Battery Hydrometer A2B Range 1.100 to 1.300 specific Gravity OM Sons / Global

Glass jar for Hydrometer

- H195 Glass jar borosilicate for Hydrometer capacity 250 ml
Certification Charges Extra As Applicable
- H196 Calibration Certificate by private for Hydrometer BAUME and Specific Gravity
- H197 Calibration Certificate charges for Brix hydro With Thermometer
- H198 Calibration Certificate charges for Brix hydrrometer / Twaddle Hydrometer
- H199 Calibration Certificate charges for Sikes hydrrometer

LACTOMETER / URINOMETER GLASS HYDROMETERS

- H200 Lactometer Tel Tru
- H201 Lactometer Liemco/Leimco Indian make / Star 0-40 / Zeal type Indian
- H202 Lactometer Star ISI RANGE 1.020-1.035/0.5 - Zeal type
- H203 Lactometer Zeal UK MAKE
- H204 Urinometer tel tru make range 1.000 to 1.060

HYGROMETERS DIAL TYPE ANALOG

- H205 Barigo German 85 mm 0 to 100% R.H. CAT NO 409
- H206 Barigo German 100 mm 0 to 100% R.H. HAIR HYGROMETER CAT NO 410
- H207 Huger 100 mm 0-100% model 3125 assemble. in India

HEAMOMETER / HEMOGLOBIN METER

Heamometer is equipped with two non-fading glass colour comparator. It is used for measuring Haemoglobin concentration in blood by Sahilis method.

- H208 Heamoglobinometer digital model 1315 Deep vision Direct Reading in Hb unit
- H209 Superior German Shali's type with round tube (Hemoglobinometer)
- H210 Top Tech Heamometer round tube, Indian
- H211 Top Tech Heamometer Plano-parallel with SQUARE tube Indian
- H212 Indian make complete set round tube model Ordinary Quality

HEMOGLOBIN ACCESSORIES

- H213 Hemoglobin Accessories Round Tube GERMAN SUPERIOR
- H214 Hemoglobin Accessories Round Tube INDIAN TOP TECH
- H215 Hemoglobin Accessories Round Tube JAYNA ORD
- H216 Hemoglobin Accessories Square Tube INDIAN TOP TECH
- H217 Hemoglobin Accessories HB Pipette GERMAN SUPERIOR
- H218 Hemoglobin Accessories HB Pipette INDIAN TOP TECH
- H219 Hemoglobin Accessories HB Pipette JAYNA ORD
- H220 Hemoglobin Accessories Rubber Tubing C Mouthpiece, Indian
- H221 Glass Stirrer (OR) Dropper (OR) Brush for Heamometer Indian Make

HEAMOCYTOMETER

Counting chamber serve to determine the number of particles per volume of blood. The Particles are

usually counted under microscope eg. Leucouyes, erthrocytes, thrombocytes etc.

- H222 Brightline Chamber Indian with RBC / WBC
- H223 Weber Bright line with German RBC/WBC Pipette
- H224 HBG Bright line chamber & German RBC / WBC Pipette
- H225 Rohem Silverlite Chamber, German RBC/WBC Pipette
- H226 Rohem Silverlite Chamber with China RBC/WBC Pipette
- H227 Rohem Silverlite Chamber with India RBC/WBc Pipette
- H228 ROHEM Roshanthal Chamber with Germen RBC/WBC Pipette

H229 China Chamber and China RBC/WBC Pipette

H230 China chamber with Indian RBC & WBC Pipette

SEDGEWICK RAFTER COUNTING CHAMBERS

H231 Sedgewick Rafter Counting Chamber (Plastic) Imported

H232 Sedgewick Rafter Counting Chamber (Plastic) Indian

HEAMOCYTOMETER SPARES : COUNTING CHAMBERS

H233 ROHEM Sperm Counting Chamber bright line Indian

H234 Roshanthal Counting Chamber bright line Indian ROHEM

H235 HBG Silverlite Counting Chamber

H236 Rohem Silverlite Counting Chamber

H237 Rohem Silverlite Counting Chamber (MOQ 50NOS)

H238 Webber silver quoted counting chamber

H239 China Bright line Counting Chamber

H240 China Neubar Counting Chamber Plain / Accutek

HEAMOCYTOMETER SPARES : RBC AND WBC PIPETTE

H241 RBC (OR) WBC pipette German price per each

H242 RBC (OR) WBC pipette China price per each

H243 Top Tech / Borosilicate Glass RBC (or) WBC Indian price per each

H244 RBC or WBC Indian Regular price per each Jayana

H245 Empty box Indian for Haemocytometer

H246 Haemocytometer cover glass 10 pair per pkt

HUMIDIFIER / FUMIGATOR

H247 Aerosal SS body 5 Litre Capacity

H248 Aerosol SS Body 5 Litre capacity with TIMER

H249 Fogtech ULV Fogger SS 304 Grade With Ss Drum 6 Ltr Digital Timer 1-99 Min

H250 Fogtech ULV Fogger Ss 304 Grade With Ss Drum 6 Ltr Mechanical Timer 1-99 Min

H251 Fogtech ULV Fogger Ss 304 Grade With Plastic Drum 6 Ltr Mechanical Timer 1-99 Min

H252 ROTASTAND FOR THE ABOVE FOGTECH FOGGER 360 DEGREE ROTATION

HUMIDITY CHAMBER REMI

H253 Model CHM-6 plus capacity 200 Temp Range 10c to 60c humidity 40% to 95% microp. Lcd version

H254 Model CHM-6 plus capacity 200 Temp Range 10c to 60c humidity 40% to 95% microp. TFT version

H255 AVAILABLE CAPACITY 280L,370L,450L,540L,800L,1000L

STABILITY CHAMBER REMI

H256 Model SC-6 plus capacity 200 Temp.10c to 60c Humidity 40% to 90% microp. LCD Version

H257 AVAILABLE CAPACITY 280L,370L,450L,540L,800L,1000L