

## **INCUBATORS DEEP VISION**

rotary ON/OFF switch, to work on 220volts AC. Temp. range 5°C above room temp. to 80°C.

Size Of Inner Chamber 300x300x300mm 12"x12"x12" Capacity 28 Ltr, 250 w Rating, 1 Shelves Analog Regulator Temp.

12"x12"x12" with Digital Temperature Controller

12"x12"x12" with Digital Temperature Controller and Fan for Circulating air

Size Of Inner Chamber 355x355x355mm 14"x14"x14" Capacity 45 Ltr, 250 w Rating, 2 Shelves Analog Regulator Temp.

14"x14"x14" with Digital Temperature Controller

14"x14"x14" with Digital Temperature Controller and Fan for Circulating air

Size Of Inner Chamber 455x455x455mm 18"x18"x18" Capacity 95 Ltr, 400 w Rating, 2 Shelves Analog Regulator Temp.

18"x18"x18" with Digital Temperature Controller

18"x18"x18" with Digital Temperature Controller and Fan for Circulating air

Size Of Inner Chamber 605x605x605mm 24"x24"x24" Capacity 224 Ltr, 750w Rating, 3 Shelves Analog Regulator Temp.

24"x24"x24" with Digital Temperature Controller

24"x24"x24" with Digital Temperature Controller and Fan for Circulating air

Optional Accessories PID Digital Temp. Controller

## **INCUBATOR REMI**

RHI-50 Hot air Incubator with 50 Liters capacity & Temperature setting range from RT + 5°C up to +65°C

RHI-80 Hot air Incubator with 80 Liters capacity & Temperature setting range from RT + 5°C up to +65°C

## **ION METER & ION ELECTRODES**

ION Analyser PH/ION Meter PH2100 Eutech (Without any Electrode)

ION Analyser PH/ION Meter ION700 Eutech (Without any Electrode)

"ION analyser ph/Ion meter ION2700 with ph electrode & atc probe but without ION electrodes RS232 port

"

ION Selective Electrodes ANTECH USA make with 500 ml Standard solution and Ion strength adjuster Price for One

(FOR FLOURIDE, NITRATE, CHLORIDE, SODIUM ETC ANY ONE)

Model 1016 Deep vision uP Based pH/ION/ MV Meter with 5 point calibration range - 1 to 14 ph  
Accuracy

± 0.0012 pH, Storage 999 ION range 0.001 to 19999 ppm Accuracy± 0.5% ( ION Electrodes not Included)

ION Meter Pen type For Sodium / Potassium / Nirtate / Calcium Anyone

### **COOLING INCUBATOR REMI PLUS SERIES**

CI-6 Plus Capacity 200 litters, Temperature range 5°C to 60°C Microprocessor ( VS-02) LED

CI-6 Plus Capacity 200 litters, Temperature range 5°C to 60°C Microprocessor ( VS-02) LCD

CI-10 Plus Capacity 280 litters, Temperature range 5°C to 60°C Microprocessor ( VS-02 ) LED

CI-10 Plus Capacity 280 litters, Temperature range 5°C to 60°C Microprocessor ( VS-02 ) LCD

CI-12 Plus Capacity 370 liters, Temperature range 5°C to 60°C Microprocessor (VS-03) LED

CI-12 Plus Capacity 370 liters, Temperature range 5°C to 60°C Microprocessor (VS-03) LCD

VS-03 Servo controlled Voltage Stabilizer

VS-02 Servo controlled Voltage Stabilizer

### **IRON CRUCIBLE**

IRON CRUCIBLE 50ML

IRON CRUCIBLE 100ML