

REFRACTOMETERS FOR SUGAR HAND HELD ERMA / DURALAB

- R1 Model A Sugar range 0-32 % for fruit, juices and soft drinks and coolant
- R2 Model B Sugar range 28-62% for sugar solution, infusions, concentrated juice
- R3 Model C Sugar range 58-90% for malt, honey, etc. extended range
- R4 Sugar Range 0-50% for sauce, ketchup, low calorie jam etc.
- R5 Sugar Range 45-82% for condensed milk, liquid sugar etc.
- R6 Range 0-18% Pocket Model or 0-10%
- R7 Range 0-80%
- R8 RANGE 0-90 %
- R9 Refractometer for Alcohol 0 - 25% / BRIX 0 - 40%
- R10 Refractometer for Battery 1.1 to 1.4

ATAGO, JAPAN MASTER M SERIES (LOW COST)

- R11A Hand Held Refractometer Range 0-33% Atago Japan Model Master M
- R11 Hand Held Refractometer Range 0-53% Atago Japan Model Master 53M
- R12 Hand held Refractometer range 28-62% Atago, Japan model Master 2M
- R13 Hand held Refractometer range 58-90% Atago, Japan model Master 3M
ATAGO JAPAN MASTER SERIES ATC MASTER T
ATC MASTER T
- R14 Hand held Refractometer range 0-33% Atago, Japan Master T
- R15 Hand held Refractometer range 28-62% Atago, Japan Master 2T
- R16 Hand held Refractometer range 58-90% Atago, Japan Master 3T
WP/ATC MASTER ALPHA
- R17 Hand held Refractometer range 0-33% Atago, Japan Master. Alpha
- R18 Hand held Refractometer range 28-62% Atago, Japan Master.2 Alpha
- R19 Hand held Refractometer range 58-90% Atago, Japan Master.3 Alpha
- R20 Hand Held Refractometer Range 0 -20% Brix Atago Japan Master 20 Alpha
MASTER H SERIES (Water Resistant, ATC, Heat Resistant), ATAGO, JAPAN
Refractometer designed for high temp samples during cooking or boiling
- R21 Model 93H range 45-93%

R22 Model 100H range 60-100%

REFRACTOMETERS DIGITAL ATAGO PAL SERIES

R23 PAL 1 range 0-53% (Resolution 0.1%) 10-75c ATC Japan Make

MERAPAL range 0-53% Made In India

R24 PAL 2 range 45-93% (Resolution 0.1%) 10-75c ATC

R25 PAL 3 range 0-93% (Resolution 0.1%) 10-75c ATC

R26 PAL 15 S for Glucose range 0 - 50 %

R27 PAL Alpha 0-85% (Resolution 0.1%) 10-75c ATC

R28 PAL Patisier BRIX & BAUME Scales (BRIX 0-85% & BAUME 0-45)

R29 PEN PRO Pen type 0.0 - 85.0 % BRIX Dip Model

R30 PAL RI Digital Refracto Meter 1.3306-1.5284 CAT NO 3850 resolution 0.0001 Accuracy: ± 0.0003 (water at 20°C)

R31 PAL BRIX / RI Range 0 - 93 % and RI 1.330 - 1.5284 CAT NO 3852

R32 PAL S Brix Range 0 - 93% Resolution 0.2% , Accuracy 0.4%

MERA PAL UREA 4695 RANGE 0-55%

REFRACTOMETER MILWAUKEE DIGITAL / HM DIGITAL

R33a Model No MA871 Brix Digital Refractometer Milwaukee Range 0-85% BRIX

R33b Model HM 0- 53 Brix Refractometer Digital Model SCM 1000

REFRACTOMETER : SALINITY

R34 Salt Range 0-100% ERMA OR DURALAB

R35 Salt Range 0 - 28% ERMA OR DURALAB

R36 Model ATC-S/Mill - Alpha Salinity Refractometer - range 0-100% with ATC ATAGO JAPAN MAKE cat no-2491

R36A Model S/Mill/M Salinity Refractometer range 0-100% non ATC Atago Japan make cat no.2493

R37 Model PAL SALT CAT NO-4250 Salinity Refractometer Digital Salt Range 0-10% (g/100g salt)- ATAGO JAPAN MAKE

Model PAL 06 SALINITY METER DIGITAL RANGE 0-100% SALINITY OF SEA WATER

R38 PAL 91 S for Ethylene Glycol (0-90%) cat no-4491

R39 Milwaukee Model MA 887 Range 0-150 ppt/ Specific Gravity 1.000 TO 1.114

REFRACTOMETER DIGITAL ABBE IMPORTED

R40 Model DR-A1 range 1,3000 to 1.7000 with built in light Digital Atago Cat no 1310

R41 Model R-9 Digital Abbe Refractometer Advance Make

ABBE REFRACTOMETER DIGITAL CHINA DURALAB REFRACTOMETER ABBE INDIAN

R42 Model R-8 Advanced make with digital temperature indicator with compact built in light system. supplied complete

with test piece, contact liquid, screw driver, dust cover Range RI Scale 1,3000 to 1.7100 accuracy 0.001/0.0005

(direct reading)& 0.001 (by estimation)

R44 Model R-4 Advance R.I 1.300 - 1.700 Accuracy 0.001 (direct) & 0.0001 (By Estimate) Sugar Scale 0 - 95%

Accuracy 1 % With Digital Temp.Indicator.

R46 Model DVAR-1 ABBE Refractometer Deep vision regular model

R47 Model 2WAJ-1 Abbe Refractometer range scale 1.3000 - 1.7000 accuracy +/- 0.0002 Digital Thermometer

screw driver, dust cover contact liquid, CHINA MAKE

REFRACTOMETER BUTYRO FOR OIL & SUGAR

R48 Advance make model R-7 for oil 0-100% and sugar 50-80% with temp, compensator

R49 ATAGO Digital Butyro Refractometer PAL BUTYRO Range 30-90% Japan make Pocket model cat no.

R49 ATAGO Digital Butyro Refractometer PR-BUTYRO Range 30-90% Japan make Portable model cat no- 3454

RAIN GAUGE

R50 Brass/Copper/SS body 5" dia and 12" height with rain collection jar and measuring jar

R51 Fibber Glass with 5" dia and 24" height with 200mm jar as per IS 5225-1969 Deep Vision make

R52 Self Recording Rain gauge with Quartz clock for 24x7 chart recording as per IS 5235-1967 for 10mm rainfall

R53 Spare measuring cylinder 200mm Borosilicate Glass

R54 Digital Electronic Cable free Rain Gauge. TFA make. Monsun

ROTARY SHAKERS / RECIPROCAL SHAKER

R55 See Orbital Shaker in INDEX "O"Page

XXX RECIPROCAL SHAKER DEEPCISION INDIAN MAKE 36"X18"X4"

RINGER TIMER

- R56 Ringer Timer MIN / SEC Big Display china
- R57 Ringer Timer Pacer 103 Loud Star Model HR
- R58 Ringer Timer Pacer 107 Big Display / Stop clock

REFRIGERATORS LABORATORY REMI

Revised On 15.06.2021

- R59 Model RLR 200 capacity 200 liters Temp. 1c to 10c electronic controller LED version
- R60 Model RLR 300 capacity 300 liters Temp. 1c to 10c electronic controller LED version
- R61 Model RLR 400 capacity 400 liters Temp. 1c to 10c electronic controller LED version

ROTARY EVAPORATOR

- R62 Rotary Evaporator with glass parts full set without Vacuum Pump DURALAB BASIC MODEL
Horizontal condensor
- R63 Rotary Evaporator with glass parts full set without Vacuum Pump DURALAB PT-65 DIGITAL RPM
motorized up/down motion

Rotary Evaporator with glass parts full set without Vacuum Pump DURALAB BASIC MODEL
Vertical condensor

Rotary Evaporator with glass parts full set without Vacuum Pump DURALAB DIGITAL RPM
Motorized Vertical condensor

Vacuum pump for above TID75S

Chiller for the above DURALAB