

STOP WATCHES

- S1 Stop Watch Digital 1/ 100 th Ordinary Low cost WINNER / KADIO & others
- S2 Stop Watch 1/100th of second digital PULSER / WINNER / ect
- S3 Stop Watch 1/100th of second Digital with date, time and one lap function with Lithum Bt. RACER Model 603 New
- S4 Stop Watch 1/100th of second Digital same as above but with 120 lap memory function - RACER Model 909
- S5 Stop Watch 1/10th of second metal body Diamond China make / Rocar

STOP CLOCKS

- S6 Stop Clock Digital New Model TOPPER T-Z
- S7 Stop Clock Digital table Model WINNER make 3in1 1/100th stop clock/count up and down timer big digit Model-107
- S8 Stop Clock Analog key winding Model indian make (SPEEDO/ VRISHI)

SPECTROPHOTOMETER VISIBLE RANGE EI-DEEP VISION

- S9 Model 301 Deep Vision Digital range 340 to 960 nm to measure % Trans. O.D. absorbent concentration
WL resolution 5mm with 3 ½ digit LED display,
- S10 Model 302 Deep Vision Digital Spectrophotometer range 340 to 960 for holding 2 square cuvette at a time provided
with a set of 2 optical glass cuvette
- S11 Model 304 Deep Vision Digital range 340-960nm Automatic CON.% Trans.Absorbance measurement 3 ½ digit data
display LED. WL in 1nm resolution.Provision for 10mm & 50mm pathlength cuvettes with two No.10mm&one no.50mm
- S12 Model 305 Deep Vision Digital range 340 to 960nm to measure % T. O.D.Con. Wavelength resolution 1nm
with 3 ½ digit LED display & 3 digit wave length display with one cuvette holds
- S13 Model 1305 Deep Vision uP Based WL range 340-690nm 16*2 line LCD display,100 data storage wavelength
resolution 1nm, printer attachment port
- S14 Model 1305 Deep Vision Spectrophotometer with RS232 interface & software

SPECTROPHOTOMETER VISIBLE IMPORTED SERIES

S15 HITEK VISION Model 3305 Wl range 340 - 1020 nm, Accuracy +/- 3 nm , 10/50 mm 4place Holder 5nm bandwidth,

S16 Model 2305 Deep vision VIS Wl range 325- 1000 nm, 4 nm bandwidth, Wl accuracy 2 nm, 4 digit LED display

10 mm path length Basic software included

S17 Model 2306 Deep vision VIS scanning Wl range 320- 1100 nm, 2 nm bandwidth, Wl accuracy 1 nm, 128x64 LCD

display Basic software included

S18

S19 Model AN-V600 ANTECH Singlebeam Visible with 1 cm glass cuvette (4 pcs) Wl 325-1000nm 5 nm BW 50mm Pathlength

S20

S21 Model AN-V1000N ANTECH Singlebeam Visible with 1 cm glass cuvette Wl 325-1000nm 4 nm BW 10mm Pathlength

SPECTROPHOTOMETER UV-VIS SINGLE BEAM EI DEEPPVISION

S21 Model 371 DV UV-Visble Range 190-1000nm, Four sample holder wavelength(uP adjustment) spectral bandwidth 2nm

S22 Model 1371 DV uP based UV-Visble Single Beam Range 190-1000nm. 16x2 LCD Display, 4 digit LED display for Wl, 4 spl.

holder of 10mm Cuvette, Stores 100. Band width 2 nm

SPECTROPHOTOMETER UV IMPORTED SERIES

S23 HITEK VISION Model 3371 UV - VIS Spectrophotometer Single Beam spectral band width 5 nm PC Connect Not Possible

S24 HITEK VISION Model 3375 UV - VIS Spectrophotometer Double Beam

S25 Model 2371 Deep vision UV/VIS single beam, Wl range 200- 1000 nm, 2 nm bandwidth, Wl accuracy 0.5 nm

128x64 LCD display, 100 mm path length optional, Basic S/W included

S26 Model 2373 Deep vision r UV/VIS single beam with scanning ,Wl range 190- 1100 nm, 2 nm bandwidth,

Wl accuracy 0.5 nm 128x64 LCD display, 100 mm path length optional Basic software included

S27 Model 2375 Deep vision UV/VIS DOUBLE beam with scanning ,Wl range 190- 1100 nm, 1.8 nm bandwidth, Wl accuracy

0.3 nm 320x240 Graphical LCD display, 8 Cell auto changer (50 mm & 100 mm path length optional) (software included)

S28 Model AN-UV-1500 ANTECH Single Beam UV/VIS with 1cm glass cuvette (4 pcs) WL 190-1000nm, 2nm Band width RS 232

S29 Model AN-UV-1500N ANTECH Single Beam UV/VIS with 1 cm glass cuvette (4 pcs) WL 190-1000nm, 2nm Band width

S30 Model AN-UV 2500 Antech Single Beam with scanning, WL range 190 - 1100, 1.8 nm bandwidth, WL accuracy 0.3 nm

320x240 Graphical LCD display, 10 mm 4 Position Cell Holder (software included)

S31 Model AN-UV 7000N Antech Double Beam with scanning, WL range 190 - 1100, 0.5, 1, 2, 4 nm bandwidth,

WL accuracy 0.3 nm 320x240 Graphical LCD display, 10 mm 4 Position Cell Holder (software included)

S32 Optional : Professional SOFTWARE for 2373

SPECTROPHOTOMETER ACCESSORIES FOR DEEP VISION

S33 Glass cuvette 10mm - (cap- 3.5 ml per pair)

S34 QUARTZ- 10 mm (3.5ml per pair)

S35 Spare Bulb 6v 1a Tungsten FUJI make

S36 Spare Bulb 6v 10w for DV Spectro

S37 Spare Lamp For Imported Spectrophotometer 12 v 20 W

SPECTRO CUVETTE DISPOSABLE

S38 Disposable square cuvettes for Spectrometer 1.5ml (100/Pkt.)

S39 Disposable square cuvettes for Spectrometer 3.5ml (100/Pkt.)

SPECTROPHOTOMETER CELLS AND CUVETTES

S40 Code MCG -1 45×12.5×3.5mm Length 01mm Volume 0.35ml

S41 Code MCG - 2 45×12.5×4.5mm Length 02mm Volume 0.7ml

S42 Code MCG - 3 45×12.5×7.5mm Length 05mm Volume 1.7ml

S43 Code MCG - 4 45×12.5×12.5mm Length 10mm Volume 3.5ml

S44 Code MCG - 5 45×12.5×22.5mm Length 20mm Volume 7ml

S45 Code MCG - 6 45×12.5×32.5mm Length 30mm Volume 10.5ml

S46 Code MCG - 7 45×12.5×42.5mm Length 40mm Volume 14ml

- S47 Code MCG - 8 45×12.5×52.5mm Length 50mm Volume 17.5ml
- S48 Code MCG - 9 45×12.5×102.5mm Length 100mm Volume 35ml
- S49 Code MCG-244 45×12.5×12.5mm Length 10mm Volume 0.7ml
- S50 Code MCG-254 45×12.5×12.5mm Length 10mm Volume 1.0ml
- S51 Code MCG 264 45×12.5×12.5mm Length 10mm Volume 1.5ml
- S52 Code MCG 594 45×12.5×12.5mm Length 10mm Volume 1.7ml
- S53 Code MCG 204 45×12.5×12.5mm Length 10mm Volume 3.5ml for flurimeter

QUARTZ (RATES PER PAIR)

- S54 Code MCQ - 1 45×12.5×3.5mm Length 01mm Volume 0.35ml
- S55 Code MCQ - 2 45×12.5×4.5mm Length 02mm Volume 0.7ml
- S56 Code MCQ - 3 45×12.5×7.5mm Length 05mm Volume 1.7ml
- S57 Code MCQ - 4 45×12.5×12.5mm Length 10mm Volume 3.5ml
- S58 Code MCQ - 5 45×12.5×22.5mm Length 20mm Volume 7ml
- S59 Code MCQ - 6 45×12.5×32.5mm Length 30mm Volume 10.5ml
- S60 Code MCQ - 7 45×12.5×42.5mm Length 40mm Volume 14ml
- S61 Code MCQ - 8 45×12.5×52.5mm Length 50mm Volume 17.5ml
- S62 Code MCQ - 9 45×12.5×102.5mm Length 100mm Volume 35ml
- S63 Code MCQ-24 45×12.5×12.5mm Length 10mm Volume 1.4ml
- S64 Code MCQ-244 45×12.5×12.5mm Length 10mm Volume 0.7ml
- S65 Code MCQ-254 45×12.5×12.5mm Length 10mm Volume 1.0ml
- S66 Code MCQ 264 45×12.5×12.5mm Length 10mm Volume 1.5ml
- S67 Code MCQ 564 45×12.5×12.5mm Length 10mm Volume 0.7ml
- S68 Code MCQ 204 45×12.5×12.5mm Length 10mm Volume 3.5ml for Flurimeter 4 sides clear
- S69 Code MCQ 214 45×12.5×12.5mm Length 10mm Volume 3.5ml for Flurimeter 4 sides clear

SEMI MICRO CELL WITH BLACK WALLS

- S70 Code MCQ-34 45×12.5×12.5mm Size 10mm Capacity 1.4ml
- S71 Code MCQ-444 45×12.5×12.5mm Size 10mm Capacity 0.45ml
- S72 Code MCQ 674 45×12.5×12.5 mm Size 10 mm Capacity 0.1ml /0.05 ml any one

SALINITY METER PEN TYPE EUTECH

S73 SALT TEST 11 Range 0-10.00 ppt Acc $\pm 1\%$ F.S. Pin IP67

S74 CST EXPERT Eutech Cond / TDS / Salinity Temp. Meter Refer P 66.

SALINITY PORTABLE EUTECH

S75 SALT603K+ Salinity Range 1 to 50 ppt Acc 0.1to5% Battery IP54

SALINITY BENCH MODEL

S76 Deep Vision 671E Salinity Meter 3 1/2 Digit LED Display with Salinity Cell Range 0-2ppt / 0-20ppt / 0-50ppt

FOR SALINITY REFRACTOMETER PLEASE SEE INDEX 'R' AND ATAGO SALT REFRACTOMETER IN 'R'

SOUND LEVEL METER

S77 SL 4001 LUTRON 35 to 130db (Auto Calibration Switch with manual adjustment)

S78 SL 4011 LUTRON 35 to 130db AC/DC Output, Data Hold

S79 SL 4012 Lutron Sound level with RS232 interface and max min readings

S80 SL 4013 range 35 to 130db with wire Probe and Rs232 interface LUTRON

S81a SL 4030 LUTRON 35 to 130db with Auto Calibration

S81b SL 4010 Lutron range 35 to 130 db 3 ranges

S82 SC 941 LUTRON 94db Sound Calibrator (single point)

S83 SC 942 LUTRON 94db Sound Calibrator (Double Point)

S84A SL4023 SD Lutron 30 to 130 db range with real time data logger with USB add RS 232 Options

S84B SL4035 SD Lutron 30 to 130 db range with real time data logger with USB add RS 232 Options

S85 Metravi 1332 Sound level meter with auto cal and filter range 35 to 130 db

S86 FUYI 824 China Sound Level Meter range 30 to 130 dB or Equivlent

SOLAR METER

S87 Solar Meter Matravi Model 207

S88 SOLAR Power meter TM 207 Digital Lutron

SOIL TESTING INSTRUMENTS

S89 Soil Thermometer in Glass Thermometer Glass range 0- 50 with 1/2 division-400 mm length

S90 Rapid Moisture Meter 0-50%

S91 VICAT Needle with Dashpot .

- S92 Pycnometer for Soil with Glass Jar and tapered cone Lid
- S93 Soil ph & Moisture meter cone type China make (ITEM AGAINST ORDER) NO GURANTEEE
- S94 Model PMS-714 Soil moisture meter with metal sharp edge probe Lutron,Taiwan
- S95 Model VG 200 Soil Moisture Meter with touch Sensor Digital Range
- S96 Soil survey meter SOIL - pH / Moisture / Temp. 4 in 1
FOR Soil pH Meter See INDEX "P"

STAGE AND OCULAR MICROMETER

- S97 Stage Erma China Micrometer
- S98 Ocular Erma China Micrometer
- S99 MICROMETER EYE PIECE 10X WITH OCULAR MICROMETER INDIAN
SILICA CRICIBLES GLOSIL MAKE

- S100 Crucible Height 25mm, Diar 41mm ,Capacity 15ml
- S101 Crucible Height 28mm, Dia 47mm, Capacity 25ml
- S102 Crucible Height 37mm, Dia 57mm ,Capacity 50ml
- S103 Crucible Height 45mm, Dia 67mm, Capacity 80ml
- S104 Crucible Height 50mm, Dia 71mm, Capacity 100ml
- S105 Crucible Height 51mm, Dia 67mm, Capacity 150ml
- S106 Crucible Height 61mm, Dia 94mm, Capacity 200ml
- S107 Lid Dia 41mm,for 15ml
- S108 Lid Dia 47mm,for 25ml
- S109 Lid Dia 57mm,for 50ml
- S110 Lid Dia 67mm,for 80ml
- S111 Lid Dia 82mm,for 100ml
- S112 Lid Dia 82mm,for 150ml
- S113 Lid Dia 92mm for 200 ml

SILICA BASIN WIT SPOUT

- S114 Silica Basin with Spout Height 23mm Dia 55 mm Capacity 20ml
- S115 Silica Basin with Spout Height 27mm Dia 75 mm Capacity 45ml
- S116 Silica Basin with Spout Height 33mm Dia 87 mm Capacity 70ml

- S117 Silica Basin with Spout Height 36mm Dia 94 mm Capacity 85ml
- S118 Silica Basin with Spout Height 36mm Dia 110 mm Capacity 100ml
- S119 Silica Basin with Spout Height 74mm Dia 113 mm Capacity 200ml

SILICA CAPSULES/ ASH CRUCIBLE

- S120 Silica Capsules Height 40mm Dia15mm Capacity 10 ml
- S121 Silica Capsules Height 48mm Dia15mm Capacity 20 ml
- S122 Silica Capsules Height 54mm Dia15mm Capacity 25 ml
- S123 Silica Capsules M.D 8ml (38mm)
- S124 Silica Capsules M.D 20ml(54 mm)
- S125 Silica capsules Rectanglar 52x42x18

VOLUMETRIC CRUCIBLE (VM)

- S126 Volumetric Crucible (VM) 25 X 38 mm ht/Lid
- S127 Volumetric Crucible (VM) spare Plunzer
- S128 Volumetric Crucible (VM) Spare Lid

SILICA BOATS

- S129 Silica Boats LXWXH 78X17X11
- Silica Boats LXWXH 103X20X12

TRIANGLES

- S130 2" Silica Triangles
- S131 3" Silica Triangles
- S132 4" Silica Triangles

STAND FOR VM CRUCIBLE

- S133 Stand 2 Holes
- S134 Stand 4 Holes

SPECTROMETER

- S135 Superior Quality Deep Vision Spectrometer deluxe with line adjustments 6"
- S136 Diffraction grating holder / prism holder

SPECTROMETER PRISM

- S137 Spectra Vision 25X25 BSC / CROWN (1.51 R.I)
- S138 Spectra Vision 25X25 DF (1.62 R.I)
- S139 Spectra Vision 25X25 EDF (1.65 R.I)
- S140 Spectra Vision 32X32 BSC / CROWN (1.51 R.I)
- S141 Spectra Vision 32X32 DF (1.62 R.I)
- S142 Spectra Vision 32X32 EDF (1.65 R.I)
- S143 Spectra Vision 38X38 BSC / CROWN (1.51 R.I)
- S144 Spectra Vision 38X38 DF (1.62 R.I)
- S145 Spectra Vision 38X38 EDF (1.65 R.I)

SODIUM VAPOR LAMP AND TRANSFORMER

SPIRIT LAMP

- S148 Stainless Steel
- S149 Glass Pyramid shape

SILICON GREASE METROARC / SILICON OIL

- S150 Silicon grease 50gm
- S151 Silicon grease 1kg packing
- S152 Aero pack Silicon oil 350 CST (1 kg)

SLIT ROD

- S153 Silicon carbide rod 8x8x300mm
- S154 Slit rod clamp s.s

SILICON RUBBER TUBING (Minimum 5 Meters only) prices given per meter.

- S157 Silicon Rubber Tubing ID 1 OD 3 Thickness 1
- S158 Silicon Rubber Tubing ID 1.5 OD 3.2 Thickness 0.85
- S159 Silicon Rubber Tubing ID 2 OD 3 Thickness 0.5
- S160 Silicon Rubber Tubing ID 2 OD 4 Thickness 1
- S161 Silicon Rubber Tubing ID 3 OD 5 Thickness 1

S162 Silicon Rubber Tubing ID 3 OD 6 Thickness 1.5

S163 Silicon Rubber Tubing ID 3 OD 7 Thickness 2

S164 Silicon Rubber Tubing ID 4 OD 6 Thickness 1

S165 Silicon Rubber Tubing ID 4 OD 7 Thickness 1.5

S166 Silicon Rubber Tubing ID 4 OD 8 Thickness 2

S167 Silicon Rubber Tubing ID 5 OD 7 Thickness 1

S168 Silicon Rubber Tubing ID 5 OD 8 Thickness 1.5

S169 Silicon Rubber Tubing ID 5 OD 9 Thickness 2

S170 Silicon Rubber Tubing ID 6 OD 8 Thickness 1

S171 Silicon Rubber Tubing ID 6 OD 9 Thickness 1.5

S172 Silicon Rubber Tubing ID 6 OD 10 Thickness 2

S173 Silicon Rubber Tubing ID 7 OD 12 Thickness 2.5

S174 Silicon Rubber Tubing ID 7 OD 13 Thickness 3

S175 Silicon Rubber Tubing ID 8 OD 10 Thickness 1

S176 Silicon Rubber Tubing ID 8 OD 11 Thickness 1.5

S177 Silicon Rubber Tubing ID 8 OD 12 Thickness 2

S178 Silicon Rubber Tubing ID 9 OD 11 Thickness 1

S179 Silicon Rubber Tubing ID 9 OD 12 Thickness 1.5

S180 Silicon Rubber Tubing ID 10 OD 12 Thickness 1.1

S181 Silicon Rubber Tubing ID 10 OD 14 Thickness 2

S182 Silicon Rubber Tubing ID 12 OD 13 Thickness 0.5

S183 Silicon Rubber Tubing ID 12 OD 17 Thickness 2.5

S184 Silicon Rubber Tubing ID 13 OD 19 Thickness 3

S185 Silicon Rubber Tubing ID 14 OD 18 Thickness 2

S186 Silicon Rubber Tubing ID 15 OD 18 Thickness 1.5

S187 Silicon Rubber Tubing ID 16 OD 22 Thickness 3

S188 Silicon Rubber Tubing ID 20 OD 27 Thickness 3.5

S189 Silicon Rubber Tubing ID 7 OD 9 Thickness 1

S190 Silicon Rubber Tubing ID 7 OD 10 Thickness 1.5

S191 Silicon Rubber Tubing ID 7 OD 11 Thickness 2

S192 Silicon Rubber Tubing ID 24 OD 32 Thickness 4

S193 Silicon Rubber Tubing ID 25 OD 30 Thickness 2.5

S194 Silicon Rubber Tubing ID 32 OD 43 Thickness 5.5

MICROLITRE SYRINGE GLASS U-TEK INDIAN

S195 Fixed Needle (FN) Capacity 5µl

S196 Fixed Needle (FN) Capacity 10µl

S197 Fixed Needle (FN) Capacity 25µl

S198 Fixed Needle (FN) Capacity 50µl

S199 Fixed Needle (FN) Capacity 100µl

S200 Fixed Needle (FN) Capacity 250µl

S201 Fixed Needle (FN) Capacity 500µl

S203 Removable Needle (RN) Capacity 5µl

S204 Removable Needle (RN) Capacity 10µl

S205 Removable Needle (RN) Capacity 25µl

S206 Removable Needle (RN) Capacity 50µl

S207 Removable Needle (RN) Capacity 100µl

S208 Removable Needle (RN) Capacity 250µl

S209 Removable Needle (RN) Capacity 500µl

S210 Removable Needle (RN) Capacity 1000µl

S211 Gas tight fixed Needle (FNGT) Capacity 5µl

S212 Gas tight fixed Needle (FNGT) Capacity 10µl

S213 Gas tight fixed Needle (FNGT) Capacity 25µl

S214 Gas tight fixed Needle (FNGT) Capacity 50µl

S215 Gas tight fixed Needle (FNGT) Capacity 100µl

S216 Gas tight fixed Needle (FNGT) Capacity 250µl

S217 Gas tight fixed Needle (FNGT) Capacity 500µl

S218 Gas tight fixed Needle (FNGT) Capacity 1000µl

S219 Gas tight Removable Needle (RNGT) 1170 Capacity 5µl

- S220 Gas tight Removable Needle (RNGT) 1171 Capacity 10µl
- S221 Gas tight Removable Needle (RNGT) 1172 Capacity 25µl
- S222 Gas tight Removable Needle (RNGT) 1173 Capacity 50µl
- S223 Gas tight Removable Needle (RNGT) 1174 Capacity 100µl
- S224 Gas tight Removable Needle (RNGT) 1175 Capacity 250µl
- S225 Gas tight Removable Needle (RNGT) 1176 Capacity 500µl
- S226 Gas tight Removable Needle (RNGT) 1177 Capacity 1000µl

HAMILTON LIQUID HANDLING PRODUCTS

700 Series Micro Litertm Syringes, 5 µl - 500 µl Fixed Needle Syringes Point style 2 Needle Length-51mm (Sharp Tip)

NOTE : NO RETURN / REPLACEMENT for HAMILTON , SO PLEASE ORDER WITH CORRECT CODE NUMBERS.

- S227 Part No. 87900 Model 75 N" Volume 5 µl Gauge 26 s
Part NO. 87919 Model 75 N Volume 5ul gauge 26s
- S228 Part No. 80300 Model 701 N" Volume 10 µl Gauge 26 s
- S229 Part No. 80400 Model 702 N" Volume 25 µl Gauge 22 s
- S230 Part No. 80500 Model 705 N" Volume 50 µl Gauge 22s
- S231 Part No. 80600 Model 710 N" Volume 100 µl Gauge 22s
- S232 Part No. 80700 Model 725 N" Volume 250 µl Gauge 22s
- S233 Part No. 80800 Model 750 N" Volume 500 µl Gauge 22s

HAMILTON LIQUID HANDLING PRODUCTS : Luer Tip (LT) syringes, Needle not included.

- S235 Part No. 80001 Model 1701 LT Volume 10 µl
- S236 Part No. 80201 Model 1702 LT Volume 25 µl
- S237 Part No.80901 Model 1705 LT Volume 50 µl
- S238 Part No. 81001 Model 1710 LT Volume 50 µl
- S239 Part No. 81101 Model 1725 LT Volume 250 µl
- S240 Part No. 81201 Model 1750 LT Volume 500 µl

HUB Needles

Mental HUB Needles for use TLL syringe LT or TLL Connectors

Note : Also available in point Size 5

FIXED NEEDLE SYRINGES (For Rheodyne Injectors) Point Style 3 needle length 51mm (Blunt tip)

S241	Part No. 87919 Model 75 N Volume 5 µl Gauge 26s
S242	Part No. 80383 Model 701 N Volume 10 µl Gauge 26s
S243	Part No. 80365 Model 701 NR Volume 10 µl Gauge 22s
S244	Part No. 80465 Model 702 NR Volume 25 µl Gauge 22s
S245	Part No. 80565 Model 705 NR Volume 50 µl Gauge 22s
S246	Part No. 80665 Model 710 NR Volume 100 µl Gauge 22s
S247	Part No. 80765 Model 725 NR Volume 250 µl Gauge 22s
S248	Part No. 80865 Model 750 NR Volume 500 µl Gauge 22s
S249	Part No. 90022 Model 22/51/2
S250	Part no. 90026 Model 26/51/2
S251	Part no. 91022 Model 22/51/3
S252	Part no. 91026 Model 22/51/3

KEL-F HUB (KF) Needles

REMOVABLE NEEDLE SYRINGES POINT STYLE 2 NEEDLE LENGTH 51 mm

S255	Part no. 87930 Model 75 RN Volume 5 µl Gauge 26s
S256	Part no. 80330 Model 701 RN Volume 10 µl Gauge 26s
S257	Part no. 80430 Model 702 RN Volume 25 µl Gauge 22s
S258	Part no. 80530 Model 705 RN Volume 50 µl Gauge 22s
S259	Part no. 80630 Model 710 RN Volume 100 µl Gauge 22s
S260	Part no. 80730 Model 725 RN Volume 250 µl Gauge 22s
S263	Part no. 90126 Model 26/51/2
S264	Part no. 90533 Model 22/51/3
S265	Part no. 90134 Model 26/51/3

Fixed Needle syringes, point style 2

LUER TIP SYRINGES, POINT STYLE 2 Luer Tip (LT) syringes, Needle not included.

S268	Part no. 80301 Model 701 LT Volume 10 µl
S269	Part no. 80401 Model 702 LT Volume 25 µl

- S270 Part no. 80501 Model 705 LT Volume 50 μ l
- S271 Part no. 80601 Model 710 LT Volume 100 μ l
- S272 Part no. 80701 Model 725 LT Volume 250 μ l
- S273 Part no. 80801 Model 750 LT Volume 500 μ l
NEEDLE LENGTH- 51mm
- S274 Part no. 80000 Model 1701 N Volume 10 μ l Gauge 26s
- S275 Part no. 80200 Model 1702 N Volume 25 μ l Gauge 26s
- S276 Part no. 80900 Model 1705 N Volume 50 μ l Gauge 22s
- S277 Part no. 81000 Model 1710 N Volume 100 μ l Gauge 22s
- S278 Part no. 81100 Model 1725 N Volume 250 μ l Gauge 22s
- S279 Part no. 81217 Model 1750 N Volume 500 μ l Gauge 22s
FIXED NEEDLE SYRNGES,POINT STYLE 2 NEEDLE LEN.51mm
- S280 Part no. 80030 Model 1701 RN Volume 10 μ l Gauge 26s
- S281 Part no. 80230 Model 1702 RN Volume 25 μ l Gauge 26s
- S282 Part no. 80930 Model 1705 RN Volume 50 μ l Gauge 22s
- S283 Part no. 80130 Model 1710 RN Volume 100 μ l Gauge 22s
- S284 Part no. 81130 Model 1725 RN Volume 250 μ l Gauge 22s
- S285 Part no. 81230 Model 1750 RN Volume 500 μ l Gauge 22s

Note: Also available in point Size 3

SOXHLET APPARATUS

- S286 Soxhlet Apparatus 3 Tests with 3 Mantles and 3 Regulators without glass parts each 250 ml Capacity
- S287 Soxhlet Apparatus 6 Tests with 6 Mantles and 6 Regulators without glass parts each 250 ml Capacity
- S288 Glass parts Soxhlet Apparatus Single set 250 ml Capacity
- S289 Glass parts Soxhlet Apparatus Single set 250 ml Capacity with boss clamp
Screw gauge and Dial Thickness Gauge and Vernier caliper and Refer Page " V "
- SUN SHINE RECORDER
- S290 Sun Shine Recorder as per National standard superior quality

STEVENSON SCREEN TO HOUSE METREOLOGICAL INSTRUMENTS

S291 Stevenson Screen Small

S292 Stevenson Screen Large

SLIDE WARMING TABLE

S293 Size 300 x 250 (10 x12) MS Top

S294 Size 300 x 250 (10 x12) SS Top

S295 Size 300 x 250 (10 x12) Alluminium Top

S296 Size 300 x 250 (10 x12) Digital