

VACUUM PUMP OIL FREE DIAPHRAGM TYPE

Oil free portable Vacuum Pressure Pump (Diaphragm Type)

A portable diaphragm type vacuum Pump cum oil-free air compressors manufactured from special grade aluminium or light weight

and strength. Built in Micro Suction Filter. Diaphragms are made of 2-play nylon reinforced neoprene rubber and valves

are of SS-316 with special sealed bearings thus eliminating lubrication. CRP Models are with Corrosive Resistant parts Extremely

useful as Vacuum Pump in laboratories, filtration instruments, gas charging, suction and labelling machines, etc. equally useful as

compressor in flame Photometers, spectrophotometer, nebuliser, laminating, press, plastic welding, oil-free spray painting, agitation,

etc. Noiseless performance. Maintenance Free.

V1 TID -15 Max flow (lts/min) 15 Max Vacuum (in. Hg.) 22" Max Press (Psig) 25 Motor HP 1/20

TID 15 WITH NRV+GAUGE WHITE BODY

V2 TID 25 S Max flow (lts/min) 20 Max Vacuum (in. Hg.) 27" Max Press (Psig) 35 Motor HP 1/16

V3 TID 25 P Max flow (lts/min) 25 Max Vacuum (in. Hg.) 22" Max Press (Psig) 25 Motor HP 1/16

TID 25 S WITH NRV+GAUGE WHITE BODY 25 Max Vacuum (in. Hg.) 22" Max Press (Psig) 25 Motor HP 1/16

V4 TID 45 S Max flow (lts/min) 45 Max Vacuum (in. Hg.) 22" Max Press (Psig) 35 Motor HP 1/8

V5 TID 45 CRP Max flow (lts/min) 30 Max Vacuum (in. Hg.) 22" Max Press (Psig) 12 Motor HP 1/8

V6 TID 75 S Max flow (lts/min) 45 Max Vacuum (in. Hg.) 27" Max Press (Psig) 60 Motor HP 1/4

V7 TID 75 P Max flow (lts/min) 75 Max Vacuum (in. Hg.) 22" Max Press (Psig) 45 Motor HP 1/4

V8 TID 75 CRP Max flow (lts/min) 75 Max Vacuum (in. Hg.) 22" Max Press (Psig) 45 Motor HP 1/4

V9 TID 100 S Max flow (lts/min) 100 Max Vacuum (in. Hg.) HP 1/2

V10 Gauge set up for the above Pumps

VACUUM PUMP WITH OIL CHINA MAKE

V12 Vacuum Pump 50 LPM With Oil China Make Compact VP 115 Model Value

V13 Vacuum Pump 50 LPM With Oil Indian make Compact with Gauge, moisture trap and Handle

VACUUM PUMP WITH OIL DEEPCISION INDIAN

V14 Capacity 50 LPM with 1/4 HP Motor Single Stage

- V15 Capacity 100 LPM with 1/2 HP Motor Single Stage
- V16 Capacity 200 LPM with 1/2 HP Motor Single Stage
- V17 Capacity 50 LPM with 1/4 HP Motor Double Stage
- V18 Capacity 100 LPM with 1/2 HP Motor Double Stage
- V19 Capacity 200 LPM with 1/2 HP Motor Double Stage

VIBRATION METER LUTRON TAIWAN MAKE

- V20 Model VB8200 Vibration meter with sensor velocity 200mm/s acceleration 200m/s' RMS, peak, RS282,hold, record. RS232 interface
- V21 Model VB8201HA Vibration meter with sensor velocity 200mm/s acceleration 200ms'RMS, peak,hold,record (Max.Min.)
- V22 Model VB8203 Vibration meter with RS-232 (Software Not Included)
- V23 Vibration meter spare Sensor
- V24 Model VB 8213 Vibration Meter Lutron
- V25 Model PVB 820 Pen Vibration Meter Acceleration 200 m/s², 20.00G Velocity 200mm/s, 20.00cm/s

VISCOMETER RED WOOD NO.1 (AS PER IP 70) DEEP VISION INDIAN

- V26 Electrically heated With energy regulator Upto 2000 CTS (IP Thermometer Optional Ref Page " T ")
- V27 Electrically heated with Digital temperature controller
- V28 Spare Cup assembly for Redwood No.1

VISCOMETER RED WOOD NO.2 (AS PER IS) DEEP VISION INDIAN

- V29 Electrically heated with energy regulator Above 2000 CTS (IP Thermometer Optional Ref Page " T ")
- V30 Spare Cup assembly for Redwood No.2
- V31 Electrically heated with Digital temperature controller

VISCOMETER SAYBOLT (AS PER ASTM D88) DEEP VISION INDIAN MAKE

- V32 Electrically heated with energy regulator (IP Thermometer Optional Ref Page " T ")
- V33 Essential Accessories: Spare Flask 60ml
- V34 Spare Heater

TAR VISCOMETER

TAR VISCOMETER

VISCOMETER BATH KINEMATICS

V35 K.V. Bath Model T 8038 with Toughend glass Window Chamber size : 12" × 9" × 12" (S.S.):
Accommodation 4 Viscometers

Temp controlled by (a) Digital temp.Controller cum indicator Acc 0.10C Temp Range +40 to 100
Deg C

V36 K.V. Bath Model T 8039 with toughened glass window Chamber size . 12" × 12" × 12" (S.S.)
Acc.6 Viscometers Temp

Controlled by (a) Digital Temp. Controller cum indicator Acc. 0.10C Temp Range +40 to 100 Deg
C

V40 Spare viscometer Holder

V41 Spare Viscometer Holder lock type SS Superior

VISCOMETER TAR (AS PER ASTM D88)

V42 Electrically heated with energy Regulator (IP Thermometer Optional Ref Page " T ")

V43 Electrically Heated with Digital Temperature Controller

VISCOMETER ROTATIONAL DURALAB

Viscometer Duralab Rotational Digital Model LI- V 200

VISCOMETER ROTATIONAL ATAGO JAPAN

V44 VISCO Portable Digital Viscometer Atago Japan (VISCO in a Box) Cat NO 6800

VISCOMETER ROTATIONAL FUNGILAB SPAIN DIGITAL

V45 VISCOMETER - PRO H For Tar with TR 11 Spindle With Thermo Spear

VERNIER CALIPER MITUTOYO

V46 Vernier Caliper 6"/150mm (.02mm/.001 inch) ordinary Cat No.531-119 w / fine

V47 Vernier Caliper 8"/200mm (.02mm/.001 inch) ordinary Cat No. 532 120

V48 Vernier Caliper 12"/300mm (.02mm/.001 inch) ordinary Cat No. 530 -119

V49 Vernier Caliper 18"/24"

V50 Digital Vernier Caliper 6"/150mm (0.01mm) mm/inch Cat No. 500-196-20

V51 Digital Vernier Caliper 8"/200mm (0.01mm)mm/inch Cat No. 500-197-30

V52 Digital Vernier Caliper 12"/300mm (0.01mm)mm/inch Cat No. 500-754-10

V53 Dial Caliper 6"/150mm (.02mm/.001 inch) Cat No.505-671

V54 Dial Caliper 8"/200mm Cat No. 505-672

V55 Dial Caliper 12"/300mm Cat No. 505-673

VERNIER CALIPER INSIZE

V56 Vernier Caliper 130mm (.02mm/.001 inch) ordinary Cat no 1233-130

V57 Vernier Caliper 180mm (.02mm/.001 inch) ordinary Cat No 1233-180

V58 Vernier Caliper 12"/300mm (.02mm/.001 inch) ordinary Cat No 1205-3002S

V59 Vernier Caliper 18"/24" (450mm or 600mm)

V60 Digital Vernier Caliper 6"/150mm (0.01mm) mm/inch Cat No 1112-150

V61 Digital Vernier Caliper 8"/200mm (0.01mm)mm/inch Cat No 1112-200

V62 Digital Vernier Caliper 12"/300mm (0.01mm)mm/inch Cat No 1112-300

V63 Dial Caliper 6"/150mm (.02mm/.001 inch) Cat no 1312-150A

V64 Dial Caliper 8"/200mm Cat no 1312-200A

V65 Dial Caliper 12"/300mm Cat no 1312-300A

V66 Data cable with Software Insize for Digital Vernier Caliper Code No 7302-SPC1A

VERNIER CALIPER CHINA

V67 Vernier Caliper 6"/150mm (.02mm/.001 inch) ordinary

V68 Vernier Caliper 8"/200mm (.02mm/.001 inch) ordinary

V69 Vernier Caliper 12"/300mm (.02mm/.001 inch) ordinary

V70 Vernier Caliper 18"/24"

V71 Digital Vernier Caliper 6"/150mm (0.01mm) mm/inch

V72 Digital Vernier Caliper 8"/200mm (0.01mm)mm/inch

V73 Digital Vernier Caliper 12"/300mm (0.01mm)mm/inch

V74 Dial Caliper 6"/150mm (.02mm/.001 inch)

V75 Dial Caliper 8"/200mm

V76 Dial Caliper 12"/300mm

SCREW GAUGE

V77 Screw Gauge 0-25 (0.01mm) / Micrometer MITUTOYO Cat No. 103-137-M-110 -25

V78 Screw Gauge 0-25 (0.001mm) / Micrometer MITUTOYO Cat No. 103-129

V79 Screw Gauge Electronic Digital Micrometer 0-25mm (.001mm) MITUTOYO Cat No. 293-821-30

- V80 Screw Gauge Mechanical counter Micrometer 0-25mm (.01mm) MITUTOYO
- V81 Height Gauge 300mm dia Digital MITUTOYO Cat No 570-312
- V82 Screw Gauge 0-25 (0.01mm) / Micrometer INSIZE Cat No. 3203-25A
Screw gauge electronic digital 0-25 cat no. 3108-25A INSIZE Reading on barrel also
- V83 Screw Gauge Electronic Digital Micrometer 0-25mm (.001mm) INSIZE 3109-25A
- V84 Screw Gauge 0-25 (0.01mm) / Micrometer CHINA
- V85 Screw Gauge Mechanical counter Micrometer 0-25mm (.01mm) CHINA
- V86 Data cable with Software Insize for Digital Micrometer Code No 7302-SPC2A

THICKNESS GAUGE DIAL /DIGITAL TYPE

- V87 Dial Gauge 10(.01) Mitutoyo 2046-S & stand 7010S MITUTOYO
Dial Gauge 10(.01) Mitutoyo 2046-S & China stand
Dial Gauge 10(.01) China with China stand
- V88 Dial thickness gauge 10mm (.01mm) MITUTOYO Cat No. 7301
- V89 Dial Thickness gauge 20mm (.01mm) MITUTOYO
- V90 Dial Thickness gauge 1mm (0.001mm) MITUTOYO Cat No. 7327
- V91 Digital Thickness gauge 12mm (0.001mm) with 10mm shaft MITUTOYO Cat No. 547-401
- V92 Digital Thickness gauge 10mm (0.01mm) with 10mm shaft MITUTOYO CAT NO 547-301
- V94 Dial thickness gauge 10mm (.01mm) INDIAN
- V95 Dial Thickness gauge 1mm (0.001mm)INDIAN
- V97 Dial thickness gauge 10mm (.01mm) INSIZE
- V98 Dial Thickness gauge 20mm (.01mm) INSIZE
- V99 Dial Thickness gauge 1mm (0.001mm) INSIZE

ULTRASONIC THICKNESS GAUGE

- V100 Digital Thickness Gauge 0-10 x 0.001 mm LC INSIZE code 2871-101
- V101 Model UTM 9 Thickness Gauge Ultrasonic Model

SKIN CALIPER / DIAL PIPE THICKNESS

- V102 SkinCaliper 10mm Insize Code 2862-101
VACUUM OVEN DEEP VISION
- V103 Vacuum Oven 12 x 12 (24 litter) with PID Controller