

WHIRLING PSYCHROMETER

- W1 Zeal Thermometers in Indian Body range 20 to 120F
- W2 Zeal Thermometers in Indian Body range -5° to 50°C
- W3 Dimple Indian make range in Centigrade -5 to 50°C - " F" Range Not available
- W4 Spare thermometer range 20 to 120F / -5° to 50°C (Set of 2 Thermometers) ZEAL MAKE

WRIST ACTION SHAKING MACHINE

- W5 It has a heavy cast iron base with rugged construction & two sides arm for holding flask.Fitted with

a F.H.P. heavy duty motor with variable speed can hold flask from 100 to 1000ml Capacity.

Supplied with cord & plug to work on 220/230 volts AC supply to Hold 8 Flask

WATER BATH SEROLOGICAL INDIAN

- W6 Serological Water Bath Plain (250mmx250mmx150mm) With Energy Regulator Single Rack
- W7 Serological Water Bath 2 racks (25cm x 17.5x17.5cm) 500 watts With Energy Regulator
- W8 Serological Water Bath 4 racks (35cm x 20 x 17.5cm) 1000 watts With Energy Regulator
- W9 Serological Water Bath Plain (250mmx250mmx150mm) Digital
- W10 Serological Water Bath 2 racks (25cm x 17.5x17.5cm) 500 watts Digital
- W11 Serological Water Bath 4 racks (35cm x 20 x 17.5cm) 1000 watts Digital

WATER BATH SS DEEP VISION

- W12 Single wall with Thermostat SS 6 Holes
- W13 Double wall with Thermostat SS 6 Holes
- W14 Double Wall with Digital Controller 6 Holes
- W15 Single wall with Thermostat SS 12 Holes
- W16 Double wall with Thermostat SS 12 Holes
- W17 Double Wall with Digital Controller 12 Holes

WATER BATH CHILLER (COOLING CIRCULATORY BATH)

- W18a Waterbath circulatory -10 to 100°C WB 2000 V
- W18b Waterbath Chiller -10 to 100°C MADE IN INDIA SUPERFIT

WATER BATH CIRUCLATORY DEEPPVISION

- W19 Circulatory Water Bath with stirrer & Digital Temp Controller Ambient to 100°C

WATER BATH SHAKER

- W20 Water Bath Shaker with Tray for 12 Flask 100 ml Model RSB 12 REMI
EXTRA TRAY FOR THE ABOVE / CLAMPS 250ML AND 500ML EXTRA COST
- W21 Water Bath Shaker for 6 Flask Deepvision make With Digital Temp. Controller
WATER STILL SS DEEP VISION
- W22 Capacity wall Type 2 Litres(1 to 5 KW)
- W23 Capacity wall Type 4 Litres(3.0 KW)
- W24 Capacity wall Type 10 Litres (8.0KW)
- W25 Capacity wall Type 10 Litres Auto cut off model (8.0)
WATER DISTILLATION GLASS DURALAB
- W26 Water Distillation All Glass on Metal Stand with Metal Heater Single 4 Lit Capacity New Cut Off Model
- W27 Water Distillation All Glass 4 Ltr Capacity on Metal Stand with Metal Heater Double New Cut Off Model
- WAX BATH / PARAFIN WAX BATH
- Paraffin Wax Bath 1LTR
- Paraffin Bath Ss Cup 6 Tu ORD With Tap

WINTROBE TUBES & STANDS

- W28 Wintrobe Superior
- W29 Wintrobe Top-Tech
- W30 Wintrobe Ordinary
- W31 Wintrobe Stand for 6 tubes
- W32 Wintrobe tube filler Indian
- WIND VANE
- W33 For wind direction roof top mounting as per IS No.5733-1070 Indian Make Direction arrows with ball bearing rotation
- WEIGHT BOXES ORDINARY
- W34 For Labs Class 100 gms BRASS/SS M1 CLASS
- W35 For Labs Class 200 gms BRASS/SS M1 CLASS
- W36 Weight Box F1 Class 23 Pcs Set Weight 1 mg to 200 gms SS304

- W37 Weight Box F1 Class 12 Pcs Set Weight 1 mg to 500 gms SS304
- W38 NABL Certificate for F1 / 200 gms weight box W36
- W39 NABL Certificate for F1 / 500 gms weight box W37
- W40 For School/College Labs Weight Box for Chemical/ Analytical Balance 100gms/200gms LOCAL CERTIFICATE for M1 class
- W41 For School/College Labs Weight Box for Chemical/ Analytical Balance 100gms/200gms NABL CERTIFICATE for M1 class
- WATER SOIL TEST KIT / MULTIPARAMETER DEEP VISION
- SEE "PAGE - P" in (pH/COND/TDS-Multiparameter)
- WIND SPEED AND WIND DIRECTION
- W42 Wind Speed,Direction and Temperature Meter Model AM 4836C,1 Mt Cable CHINA
- W43 RS232C cable &software for above
- WEATHER STATION
- W44 Weather Station,TFA make Model NEXUS with wireless sensors
- W55 Digital Electronic cable free weather station OREGON scientific make WMR 200A
(Installation Charges For Weather station Within Tamilnadu Rs- 5000/-)
- W46 Walking Cold Cabinets REMI 5000 Liters to 50000 Liters Capacity