Bottle Top Dispenser







Fisherbrand[™] offers three lines of bottle top dispensers which caters to every laboratory need from simple dispensing to handling of aggressive and hazardous chemicals and reagents. Fisherbrand dispensers are robust in design with a PTFE piston for smooth movement, easy volume setting, high chemical compatibility and an adjustable delivery nozzle for ease of working in challenging laboratory conditions. These dispensers are fully autoclavable and have easy user re-calibration.

Bottle Top dispenser with re-circulation valve facilitates quick and definite purging (removal of air) without loosing any reagent.

Bottle Top Dispenser with Dual Inlet facilitates multiple dispensing modes including rinsing, dilution and bottle re-filling without dismounting the instrument from the reagent bottle.

Whatever your application may be, you can always select a suitable Fisherbrand Bottle Top Dispenser.

Features	Bottle Top Dispenser	Bottle Top Dispenser with re-circulation valve	Bottle Top Dispenser with dual inlet	
Dual Inlet Technology			✓	
Re-circulation Valve		√	✓	
PTFE Piston	✓	√	√	
180° Turn Volume Adjustment Knob	√	√	√	
Adjustable Delivery Nozzle	✓	√	✓	
360° Rotation of Dispenser on reagent bottle		√	✓	
Individual calibration report with each unit in accordance with ISO 8655 standards.	✓	√	√	
Five Bottle Adapters (28,32,36,40,45mm)	✓	√	✓	



Specifications & Ordering Information

Cat. No.	Vol. Range	Increment	Acc ±%	uracy ± ml	±%	CV ± ml	
Bottle Top	Dispenser						
13-615-10	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005	
13-615-11	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010	
13-615-12	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020	
13-615-13	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060	
13-615-14	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120	
13-615-15	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200	
Bottle Top Dispenser with Re-circulation Valve							
13-681-524	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005	
13-681-525	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010	
13-681-526	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020	
13-681-527	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060	
13-681-528	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120	
13-681-529	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200	
Bottle Top Dispenser with Dual Inlet							
13-681-530	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005	
13-681-531	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010	
13-681-532	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020	
13-681-533	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060	
13-681-534	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120	
13-681-535	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200	

In the United States:

For customer service, call 1-800-766-7000 To fax an order, use 1-800-926-1166 To order online: fishersci.com

In Canada:

For customer service, call 1-800-234-7437 To fax an order, use 1-800-463-2996 To order online: fishersci.ca

