

PRODUCT DESCRIPTION FOR "STERILE GARMENT STORAGE CABINET"

PRODUCT DESCRIPTION	
Cleanliness	Class 100 / ISO CLASS 5 / M 3.5 AREA
Velocity at Filter	90 ± 20 FPM
MOC	<ul style="list-style-type: none"> Main Body - SS-304 Door - Swing Type Door with accessories – 01 Nos. The Cabinet's inside bottom plate is S.S perforated tray.
Supply Filter	<ul style="list-style-type: none"> Micro Fiberglass Paper Media Mini-pleated HEPA Filter down to 0.3 µ Particulate at an efficiency of 99.997% duly tested on PAO/DOP tested. PAO/DOP test port at upstream of HEPA Filter. Make: Hygieno Filters
Exhaust Filter	<ul style="list-style-type: none"> Non-woven synthetic media Pre-Filter down to 10 µ Particulate at an efficiency of 90%. Washable type. Make: Hygieno Filters
Fresh Air Filter	<ul style="list-style-type: none"> Non-woven synthetic media Pre-Filter down to 5 µ Particulate at an efficiency of 75 to 80%. Washable type. Make: Hygieno Filters
Blower Assembly	<ul style="list-style-type: none"> Dynamically and statically balanced motor and blower.
Electrical	<ul style="list-style-type: none"> Feather Touch On/Off Panel LED Light – 1 No. (Make: Phillips or Halonix or equivalent) UV Light – 1 No. (Make: Phillips or equivalent)
Pressure Measurement	<ul style="list-style-type: none"> 1 No. Magnehelic Pressure Gauge to check the pressure drop across the HEPA filter
Power Supply	<ul style="list-style-type: none"> 230 V, 50Hz, Single Phase
Documentation	<ul style="list-style-type: none"> DQ, IQ, OQ (Our Factory Tested)

