

SPECTROPHOTOMETER CM-3700A

THE NEW LEADER

Konica Minolta's flagship model: High accuracy and repeatability even at low reflectance/transmittance



Giving Shape to Ideas

A high-accuracy, top-class model filled with Konica Minolta's advanced optical technology. Ideal for demanding customers with strict requirements for R&D, QC, and CCM applications.

An advanced spectrophotometer for reflectance and transmittance measurements of a broad range of subjects including pigments, dyes, plastic, textiles, paints, ceramics, etc.

Top-of-the-line model providing high accuracy and high reliability

Konica Minolta's advanced optical, sensing, and signalprocessing technology provide excellent repeatability.

Strictly selected highquality parts ensure long-term stability and reliability.

Strict accuracy control traceable to national standards ensures high quality with high inter-instrument agreement.

56 Polychromatic unit

Illumination/viewing geometry conforms to CIE, ISO, ASTM, DIN, and JIS requirements for di:8°, de:8° (diffuse illumination/8° viewing) geometry for reflectance and to CIE, ASTM, DIN, and JIS requirements for transmittance.

2 Pulsed xenon arc lamp light source provides high stability, long life, and excellent repeatability even on dark and highchromaticity colors.

3 6-inch integrating sphere has a powdered barium sulfate (BaSO,) coating with superior optical characteristics.

4 Double-beam feedback system directly monitors the light emitted by the xenon lamp at the time of measurement and automatically compensates for changes in brightness or spectral characteristics to ensure highaccuracy measurements.

A diffraction grating ^⑤ provides highefficiency separation of light by wavelength, ensuring high repeatability even when measuring dark colors, and the silicon photodiode array sensor [©] quickly converts the separated light into electrical currents. These elements are mounted on stainless steel with a low coefficient of thermal expansion to ensure long-term stability.

The optical systems for sample measurement and light-source monitoring are mounted on an aluminum alloy block for long-term stability.

* CM-3700A is computer-controlled. Software such as optional SpectraMagic NX required





Top standard instrument for color control systems

With its high accuracy, high repeatability, and high reliability, the CM-3700A can be used as the top standard instrument to which other instruments are referenced in color quality control systems including multiple instruments, both within a company or between companies.



Strict quality assurance system ensures reliability and peace of mind

Konica Minolta's color-measuring instruments are traceable to national standards for wavelength accuracy and white calibration for strict accuracy control. Konica Minolta has also received ISO 9001 and ISO 14001 certification for its integrated guality control system from product development through manufacturing and after-sales service.

Comprehensive measurement functions

Reflectance measurements

8 Switchable between SCI and SCE measurements

SCI (specular component included) measurements minimize the influence of surface conditions on measured values, making it suitable for CCM applications. SCE (specular component excluded) measurements correspond closely to professional visual evaluation.

9 Selectable measurement areas

Measurement areas of Ø25.4mm, Ø8mm, and 3×5mm (rectangular) can be selected according to the application.

10 Variable UV for measuring fluorescent samples

The UV cut filter can be adjusted in 1000 steps for measurements of samples containing optical whitening agents such as paper, pulp, etc.



Transmittance measurements

Ð Sideless transmittance chamber for unlimited sample length (Maximum thickness: Approx 50 mm)

di:0°, de:0° (diffuse illumination/0° viewing angle) geometry for spectral transmittance measurements of sheet samples such as glass, filters, etc. as well as cells containing liquid samples such as foods, cosmetics, etc



Additional Features of CM-3700A

Long sample holder arm enables measurement at center of A4-size sample.

open when opened fully for easy positioning of thick samples.





Comprehensive support system

Konica Minolta has service facilities worldwide to perform quick inspection, calibration, and repair and ensure that your instruments always provide their best performance.

Sample holder arm stays

Frame around power switch prevents switching power off by mistake.





USB communication interface for easy connection and highspeed communication



System Diagram



Optional Accessories

Transmittance Specimen Glass Cell Holder CM-A96 CM-A97/CM-A98/CM-A99



Calibration Plate CM-A100

SpectraMagic NX (Optional)

SpectraMagic NX enables you to perform comprehensive color inspection and analysis of incoming raw materials, in-process production, and outbound color-critical goods and materials in virtually any industry. With SpectraMagic NX you can insert digital images with measured data. Measure samples in any of 8 universally accepted color spaces. Select from 16 illuminants, and up to 40 indices to determine specific color and appearance properties, such as strength, brightness, haze, yellowness, opacity and whiteness. You can even configure up to 8 customized color equations. Reports range from simple Pass/Fail to trend charts, histograms, color plots, and spectral graphs. SpectraMagic NX comes with predefined templates, or you can create your own templates. For illustrations and explanations to understanding color and color measurement technology, there is a link to Konica Minolta's wellknown and respected "Precise Color Communication", as well as step by step navigation help. Available in 8 languages: English, French, German, Italian, Spanish, Japanese, Chinese (Simplified/Traditional), and Portuguese. * Windows® is a trademark of Microsoft Corporation in the USA and other countries.

Supports Windows[®] 7/8.1/10

Transmittance Zero



Color Plates (14 colors)



Specifications

Illumination/viewing system	Reflectance	di:8°, de:8° (diffuse illumination/8° viewing angle)			
		SCI (specular component included)/SCE (specular component excluded) switchable			
		Conforms to CIE No.15, ISO 7724/1, ASTM E 1164, DIN 5033 Teil 7 and JIS Z 8722 condition c standard.			
	Transmittance	di:0°, de:0° (diffuse illumination/0° viewing angle)			
		Conforms to CIE No.15, ASTM E 1164, DIN 5033 Teil 7 and JIS Z 8722 condition g standard.			
Detector	Silicon photod	ode array with flat holographic grating			
Wavelength range	360 to 740 nm				
Wavelength pitch	10 nm				
Half bandwidth	Approx. 14 nm average				
Measuring range	0 to 200%; Resolution: 0.001%				
Light source	Pulsed xenon arc lamp				
Integrating sphere size	ø152 mm				
Minimum measurement interval	3 seconds				
Measurement/illumination area	Reflectance:C SAV : 3x5 m MAV : ø8 mr LAV : ø25.4 Transmittance:	hangeable between SAV, MAV, and LAV Im measurement / Ø5X7 mm illumination m measurement / Ø11 mm illumination mm measurement / Ø28 mm illumination Approx. ø20 mm / Ø25 mm			
Repeatability	When white calibration plate is measured 30 times at 10-sec. intervals after white calibration has been performed: Spectral reflectance: Standard deviation within 0.05% Chromaticity: Standard deviation within ΔE^*_{ab} 0.005				
	When black tile (BCRA Series II; reflectance: 1%) is measured 30 times at 10-second intervals after white calibration has been performed: Spectral reflectance: 380 to 740 nm: Standard deviation within 0.02% 360 and 370 nm: Standard deviation within 0.04% Chromaticity: Standard deviation within ΔE^*_{ab} 0.05				
Inter-instrument agreement	Mean $\Delta E^*_{ab} 0.0$ Max $\Delta E^*_{ab} 0.3$ II color tiles co	08 (typical) Average for 12 BCRA Series II color tiles. (corresponds to approx. ΔE^*_{CMC} 0.2) for any of 12 BCRA Series mpared to values measured with Konica Minolta master body.			
UV adjustment	Computer cont	trolled: continuously variable			
Transmittance chamber	Maximum sam Maximum sam cover is open) Sample holder samples can b	ple thickness: Approx. 50 mm ple length: Unlimited (no sides when transmittance chamber (optional) for holding sheet samples or containers of liquid e installed/removed			
Interface	USB 1.1				
Power	AC 100 to 240	V 50/60 Hz 25 VA (using included AC Adapter)			
Operation temperature/ humidity range (*1)	13 to 33°C, rel	ative humidity 80% or less (at 33°C) with no condensation			
Storage temperature/ humidity range	0 to 40°C, rela	tive humidity 80% or less (at 35°C) with no condensation			
Size (WxHxD)	271 x 274 x 50	0 mm (10-11/16 x 10-3/4 x 19-11/16 in.)			
Weight	18 kg (39.7 lb.				

*1 Operation temperature/humidity range of products for North America : 13 to 33°C, relative humidity 80% or less (at 31°C) with no condensation

- Displays shown are for illustration purposes only.
 The specifications and appearance shown herein are subject to change
- without notice.
- Company names and product names used herein are trademarks or registered
- registered trademarks of their respective companies. KONICA MINOLTA, the Konica Minolta logo and symbol mark, "Giving Shape to Ideas", and SpectraMagic are registered trademarks or trademarks of KONICA MINOLTA, INC.

KONICA MINOLTA, INC. Konica Minolta Sensing Americas, Inc. Konica Minolta Sensing Europe B.V.

Konica Minolta (CHINA) Investment Ltd.

Konica Minolta Sensing Singapore Pte Ltd.

Osaka, Japan New Jersey, U.S.A. European Headquarter /BENELUX German Office French Office UK Office Italian Office Swiss Office Nordic Office Polish Office Turkish Office SE Sales Division Beijing Office Guangzhou Office Chongqing Office Qingdao Office Wuhan Office

Phone : 888-473-2656 (in USA), Nieuwegein, Netherlands München, Germany Roissy CDG, France Warrington, United Kingdom Cinisello Balsamo, Italy Dictition, Switzodond Västra Frölunda, Sweden Wroclaw, Poland Istanbul, Turkey Shanghai, China Beijing, China Guangdong, China Chongqing, China Shandong, China Hubei, China Singapore Goyang-si, Korea

201-236-4300 (outside USA)
Phone: +31(0)30 248-1193
Phone: +49(0)89 4357 156
Phone: +33(0) 1 80 11 10 7
Phone: +44(0) 1925 46730
Phone : +39 02849488.00
Phone: +41(0)43 322-9800
Phone: +46(0)31 7099464
Phone: +48(0)71 73452-1
Phone : +90 (0) 216-528 56
Phone : +86- (0)21-5489 02
Phone : +86- (0)10-8522 15
Phone : +86- (0)20-3826 42
Phone : +86- (0)23-6773 49
Phone: +86- (0)532-8079 1
Phone : +86- (0)27-8544 99
Phone : +65 6563-5533
Phone: +82(0)2-523-9726

Worldwide sales network

ISO Certifications of KONICA MINOLTA, Inc., Sakai Site				
ISO 9001	ISO 14001			
JQA-QMA15888	JQA-E-80027			
Design, development, manufacture/	Design, development,			
manufacturing management, calibration, and	manufacture, service and sales			
service of measuring instruments	of measuring instruments			

SA)	Fax :	201-78	5-24	82		
1193	Fax :	+31(0)	30 2	4 81	211	
7 156 0	Fax :	+49(0)	89 4	357	156 9	99
10 70	Fax :	+33(0)	1 80	11	10 82	
67300	Fax :	+44(0)	1925	5 71	1143	
.00	Fax :	+39 02	8494	188.	30	
9800	Fax :	+41(0)	43 3	22-9	809	
9464						
52-11	Fax :	+48 (0)717	734 :	52 10	
28 56 56	Fax :	+90(0)	212	-253	3 49 6	9
39 0202	Fax :	+86-(0)21-	5489	0005	ō
2 1551	Fax :	+86-(0)10-8	8522	2 124	1
26 4220	Fax :	+86-(0)20-3	3826	6 4223	3
'3 4988	Fax :	+86-(0)23-6	6773	3 4799	Э
079 1871	Fax :	+86-(0)532	-807	9 187	73
4 9942	Fax :	+86-(0)27-8	8544	999	1
3	Fax :	+65 65	60-9	721		
9726	Fax :	+82(0)	31-9	95 - 6	511	
	 -		-		-	

https://konicaminolta.com/instruments/network

Konica Minolta Sensing Korea Co., Ltd.

9242-4851-32 BJDPK13

Main service facilities



SAFETY PRECAUTIONS For correct use and for your safety, be sure to read the

instruction manual before using the instrument. Always connect the instrument to the specified power supply voltage. Improper connection may cause a fire or electric shock.