

The Leneta Company

15 Whitney Road

Mahwah, NJ 07430-3129 USA

Phone: 201.847.9300 • **Fax:** 201.848.8833

Email: sales@leneta.com • **Website:** www.leneta.com

Item # 5C, Form 5C Opacity Chart



Form 5C Opacity Chart

Form 5C is a black and white sealed opacity chart. The top and bottom are white and the middle section is black. Overall dimensions are 7-5/8 x 10-1/4 in (194 x 260 mm).

The header includes a writable area for recording:

Date, Time

Formula

Batch Number

Notes

A 1/4" diameter hole is centered just under the top of the chart.

SPECIFICATIONS

Color	Black & White
Size	7-5/8 x 10-1/4 in 194 x 260 mm
Box Quantity	250
Boxes per Case	4
Box Weight	9 lb