

Universal Statistics Printer



FEATURES

- **Easy CONNECT**
No Skilled Engineer required
- **User Friendly Weighing Balance**
 - User can start using in Minutes, No need to pay Extra Installation & Demo Cost
 - Universal Printer to Connect with Any Brand of Weighing Balance having built-in Printer Interface
- **Built-in Clock**
 - Enables Printing Time & Date information even with a balance which doesn't have a built-in clock function
- **Enhanced Support for ISO/GLP**
- **Protect Date & TIME Setting by Password**
 - Prevents Unintended changes of time and date setting
- **Automatically Sets the appropriate communication with a balance**
- **Brightly illuminated OLED Dot Matrix Display for Easy-to-read**
- **Standard Plain Paper Print out includes :**
 - Company Name, Balance Model, Balance Serial Number, ID, Date, Start Time, Sample No., #Sample Readings, End Time, Signature
 - Operator can choose items to Print Out
- **Various Modes Expand Utility**
 - Statistical Calculation Mode (Results #Samples, Total Samples, Max., Min., Range, Mean Value, Standard Deviation, Coefficient Variation)
 - Tare/Net/Gross Weight Printing
 - SUM Mode (Print out the total values of samples)
 - Averaging Mode
 - Formulation Mode (Print out each sample's weight independently when blending samples)
 - Comparator Mode (Print out results whether weight of sample is within target range)
 - Pipette Calibration (Calibrate Capacity of Pipette & Flask)
 - Interval Printing Mode (Print out current weight at a regular interval)
- **Affordably Priced**
Economically Priced Consumables Readily Available. e.g. Paper Roll, Ink Ribbon

ORDERING INFORMATION

Order Code	Description
USP100	Universal Statistics Printer