



Personalized PDF Catalog
Catalogue Generated 11 октября 2023 г.

Grade 0858 seed test paper, Whatman™

Supplier: Whatman products (Cytiva)






Grained, with medium fast flow rate and medium retention. A universal filter paper, used for the filtration of extracts, oils, beer, syrups etc.; also applied in filter presses or for the aspiration of liquids.



These filter papers are made from super-refined cellulose and have been specifically designed to have particular properties for each technical application. The applications range from the filtration of beverages to the purification of electroplating baths.

Available prepleated as Grade 0858 ½.

ORDER

VWR Catalog Number	Quantity	Grade	Size	Pk
515-8015		0858	110×580 mm	500
515-8001		0858	390×390 mm	500
515-8002		0858	450×450 mm	500
WHAT10334346		0858	570×580 mm	100
WHAT10390046		0858	1060×560 mm	100

SPECIFICATIONS

Grade	0858
Weight	75 g/m ²
Thickness	170 µm
Particle retention	7 - 12 µm
Filtration speed	55 s (Herzberg)